
CE

CINCH<sup>®</sup> 930LCINCH<sup>®</sup> 417M

## Tube Securement Device

ES	Dispositivo de aseguramiento de sonda
DE	Schlauchsicherungseinheit
IT	Dispositivo di fissaggio del tubo
FR	Dispositif de fixation de tube
NO	Slangefesteeinheit
SV	Fästanordning för sond
PT-BR	Dispositivo de Fixação do Tubo
NL	Buisbevestigingsapparaat
DA	Slangesikringsenhed
FI	Putken kiinnityslaite
AR	جهاز تثبيت الأنبوب
BG	Устройство за захващане на тръба
CS	Zařízení k zajištění sondy
ET	Toru kinnitamise seade
EL	Συσκευή ασφάλισης σωλήνα
HU	Szondarögzítő eszköz
JA	チューブ固定デバイス
LV	Caurules nostiprināšanas ierīce
LT	Vamzdelio sutvirtinimo įtaisas
PL	Urządzenie do mocowania zgłębnika
PT	Dispositivo de Fixação do Tubo
RO	Dispozitiv de asigurare tub
SK	Zariadenie na zabezpečenie trubičky
SL	Pripomoček za pritrditev cevke
TR	Tüp Sabitleme Cihazı

C4732-A 12/2020

Made in U.S.A.

HYPOALLERGENIC

DEHP

Tube Securement Device

Instructions for Use

**Intended Use**

The tube securement device is intended to provide external securement of medical tubes and lines.

**Indications for Use**

The tube securement device is indicated for use in patients of all ages who require external securement and strain relief of medical tubing.

	Tubing size range
CINCH® 930L	9F - 30F
CINCH® 417M	4F - 17F

**Contraindications**

Contraindications for the tube securement device include but are not limited to:

- Known allergies to silicone or silicone based adhesive
- Compromised skin conditions such as open wounds / sores, rash, infection, or burns
- Placement closer than 5 cm to any tubing insertion site

**Clinical Benefits**

Clinical benefits to be expected when using the tube securement device include but are not limited to:

- Prevents lines from becoming tangled or snagged
- Provides motion control where the lines enter the body
- Provides external strain relief of medical tubing
- Initial defense for inadvertent pulls that could otherwise displace the line

**Performance Characteristics**

Performance characteristics of the tube securement device include but are not limited to:

- Accommodate pertinent tubing size
- Comfortable when not in use
- Permit multiple external securement cycles
- Non-irritating to skin for majority of patient population

**Placing the tube securement device**

**Note:** This device is intended to be used by qualified professionals and lay users.

**Note:** Remove tube securement device from backing. Tube securement device should be placed on clean, healthy, intact skin at a minimum of 5 cm from any tube insertion points.

Clean Surface

Apply device to surface

Place tube on device (perpendicular to the strap)

Pull the tip of the strap through the slit in the device  
**Caution: Do not overtighten**

**Caution:** If redness or irritation occurs, stop use and consult a healthcare professional.

**Device Longevity**

The tube securement device is meant to be periodically replaced for optimal performance, functionality, and cleanliness. Typical device longevity ranges from 3-5 days. The tube securement device should be replaced if tears, adhesive is separating from the skin, or other signs of failure are noted. To replace the device, first remove the device by releasing the clasp mechanism and gently removing the device from the skin. The device can be placed in conventional waste for disposal. The skin should be cleaned and a new device applied if needed.

**WARNING:** This device is intended for single use. Do not reuse, reprocess or re-sterilize this medical device. Doing so may compromise biocompatibility characteristics, device performance and/or material integrity; any of which may result in potential patient injury, illness and/or death.

**Note:** Please contact AMT, our Authorized Representative (EC Rep), and/or the competent authority of the member state in which you are established if a serious incident has occurred in relation to the device.

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Manufactured for Applied Medical Technology, Inc.