# AMT Traditional G-JET®

Traditional Length GJ Feeding Device

## Directions for Use - Tube Use and Care Instructions



For placement instructions for hospital professionals please see [Device Placement Directions for Use](#) document.

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Made in the U.S.A.  
C4688-E 07/2024

Read the following warnings and instructions for use prior to using the Traditional G-JET® device. Not following these guidelines can result in device failure and/or serious injury and/or death. The Traditional G-JET® is considered safe for the intended use when used within the guidelines listed throughout these instructions.

**PRIOR TO USE:**

- Components are supplied sterile: for single user only. Do not re-use or re-sterilize.
- Prior to placement, please inspect all contents of the kit for damage. If the package is damaged or sterile barrier is breached, do not use the product.

**DEVICE PLACEMENT AND REPLACEMENT:**

- Due to the presence of the tubing support (spring) in the tubing, do not cut the device to custom length. Do not cut the distal end of the tubing to create a customized jejunal length. Doing so will eliminate the soft, tapered tip of the device and expose the tubing support.
- The suture loop at the distal end of the device is most often utilized during endoscopic placements. If the suture loop of the device is not necessary for the anticipated placement method, it is recommended to cut and gently remove the suture loop prior to placement. Do not use excessive force to remove the suture loop as this can damage the tapered end of the tubing.
- A gastropexy must be performed to affix the stomach to the anterior abdominal wall, the feeding tube insertion site identified, stoma tract dilated and measured prior to initial tube insertion to ensure patient safety and comfort. The length of the tube should be sufficient to be placed beyond the ligament of treitz. Do not use the retention balloon of the feeding tube as a gastropexy device. The balloon may burst and fail to attach the stomach to the anterior abdominal wall.
- Removal is recommended to be performed by a qualified clinician. Spontaneous closure of stoma may occur as early as twenty-four (24) hours after removal. Insert a new device if enteral feeding by this route is still intended. If closure is desired, apply a dressing over the stoma site.

**GENERAL USE:**

- Do not exceed the maximum rated balloon volume. Do not use air. Do not inject contrast into the balloon.
- Never connect the jejunal port to suction. Do not measure residuals from the jejunal port. Do not use continuous or high intermittent suction. High pressure could collapse the tube or injure the stomach tissue and cause bleeding.
- Do not channel medications through the jejunal feeding port. This will clog the tube. To help avoid clogging, use a feeding pump to deliver formula through the jejunal lumen. Never allow formula to stand in the tube.
- This device is intended to connect to compatible enteral devices only. Do not use for non-enteral applications.
- Do not use excessive force or pressure when channeling feed or medication, or to attempt to clear a clog in the tubing. This can cause the tubing to rupture or lead to damage of the tubing support structure. If clog cannot be cleared or tubing is clogging on a frequent basis, this may indicate that the device needs replaced. Failure to follow this warning may result in device failure and/or serious patient injury.
- Do not use excessive force to flush the tube. Excessive force can perforate the tube and can cause injury to the gastrointestinal tract. **Use a 30 to 60 ml syringe. Do not use smaller size syringes as this can increase pressure on the tube and potentially rupture smaller tubes.**
- When using a bolus style connector, this device has the potential to misconnect to intravenous, breathing system, limb cuff, and temperature connectors.
- This device is intended for single use. Do not reuse, re-sterilize or reprocess this medical device. Doing so may compromise biocompatibility characteristics, device performance and/or material integrity; any of which may result in potential patient injury, illness and/or death.

**WHEN USING A TRANSITION ADAPTER:**

- This device has the potential to misconnect with small bore connectors of other healthcare applications. Use this device only to connect to compatible enteral devices. Do not use for non-enteral applications.
- Tighten only by hand. Never use excessive force or a tool to tighten a rotating connector. Improper use can lead to cracking, leakage, or other failure.
- Ensure device is connected to an enteral port only and NOT to an IV set.

**BALLOON INFLATION VOLUME**

Inflate the balloon with sterile or distilled water within the fill-volume range below:

Table 1			
FR Size	Min Volume	Recommended Volume	Max Volume
14 FR & 16 FR	3 ML	4 ML	5 ML
18 FR	7 ML	8 ML	9 ML
22 FR	8 ML	10 ML	12 ML

- 14Fr and 16Fr sized tubes: Inflate the balloon with 3-5 ml of sterile or distilled water.  
**CAUTION: Do not exceed 5 ml total balloon volume. Do not use air. Do not inject contrast into the balloon.**
- 18Fr sized tubes: Inflate the balloon with 7-9 ml of sterile or distilled water.  
**CAUTION: Do not exceed 9 ml total balloon volume. Do not use air. Do not inject contrast into the balloon.**
- 22Fr sized tubes: Inflate the balloon with 8-12 ml of sterile or distilled water.  
**CAUTION: Do not exceed 12 ml total balloon volume. Do not use air. Do not inject contrast into the balloon.**

**KIT CONTENTS**

**CAUTION: Components are supplied Sterile: For single user only. Do not re-use or re-sterilize. For enteral nutrition and /or medication only.**

KIT CONTENTS	
(1) Traditional Length GJ feeding device	(1) Small Syringe (Luer Slip)
(1) Large Syringe (Catheter or ENFit®)	(2) 4 "x 4" (10 cm x 10 cm) Gauze

**DESCRIPTION**

The AMT Traditional Length GJ feeding device (**Fig. A**) provides for simultaneous gastric decompression / drainage and delivery of enteral nutrition into the distal duodenum or proximal jejunum. It enters the stomach through a gastric stoma. The tube is held in place (within the stoma tract) by means of an inflatable balloon (**Fig. A-1**) and a sliding external bolster (**Fig. A-2**). The external tri-port contains two large ports; one labeled "JEJUNAL" and one labeled "GASTRIC". The JEJUNAL port is used for feeding into the small intestine. The GASTRIC port is used to drain the stomach with the use of low intermittent suction or gravity drainage. A third port labeled "BAL" is used to inflate and deflate the balloon.

**INTENDED USE**

The AMT Traditional G-JET<sup>®</sup> is intended for gastric decompression / drainage and delivery of enteral nutrition into the distal duodenum or proximal jejunum. The AMT Traditional G-JET<sup>®</sup> is intended to be placed by qualified healthcare professionals. The AMT Traditional G-JET<sup>®</sup> is intended to be used by clinicians and caregivers/users.

**INDICATIONS FOR USE**

The Traditional Length GJ Feeding Device is indicated for use in adult, adolescent, child, and infant patients over 10kg who cannot absorb adequate nutrition through the stomach, who have intestinal motility problems, gastric outlet obstruction, severe gastroesophageal reflux, are at risk of aspiration, or in those who have had previous esophagectomy or gastrectomy. The use of this tube is also clinically indicated when simultaneous gastric decompression and jejunal feeding are needed. This includes patients in whom malnutrition already exists, or may result, secondary to concurrent conditions.

**CONTRAINDICATIONS**

Contraindications for placement of a transgastric-jejunal feeding device include, but are not limited to ascites, colonic interposition, portal hypertension, peritonitis and morbid obesity.

**COMPLICATIONS**

The following complications may be associated with any transgastric-jejunal feeding device: Skin Breakdown • Infection • Hypergranulation Tissue • Stomach or Duodenal Ulcers • Intraperitoneal Leakage • Pressure Necrosis • Jejunal Limb Migration • Bowel Perforation • Tubing Support (Spring) Dislodgement • Intussusception

**NOTE:** Verify package integrity. Do not use if package is damaged or sterile barrier compromised.

**NOTE:** Please contact AMT, our Authorized Representative (EC Rep), and/or the competent authority of the member state in which you are established if a serious incident has occurred in relation to the device.

**CLINICAL BENEFITS**

Clinical benefits to be expected when using the AMT Traditional G-JET<sup>®</sup> include but are not limited to:

- Provides direct access to the stomach and jejunum for feeding, decompression, or channeling medication
- Device may be placed either during an initial placement procedure or as a replacement device
- Channeling of medication and feeding to the patient is immediate and sustains life
- Made from medical-grade silicone to reduce skin irritation and improve patient comfort
- Adjustable bolster to accommodate various patient sizing needs
- Allows patients to feed longer through less instances of the tubing kinking or clogging

**PERFORMANCE CHARACTERISTICS**

Performance Characteristics of the AMT Traditional G-JET<sup>®</sup> include but are not limited to:

- Unique internal tube structure – transitions from a tri-lumen design in its proximal (gastric) segment to one functional lumen in its distal (jejunal) segment
- Jejunal lumen offers a larger cross-sectional area in the jejunal portion to help prevent clogging
- Anti-kink technology (16F and 18F only) to help maintain flow
- Unique "apple" shape balloon designed to reduce leakage and reduce the risk of accidental pull-outs from the stoma site
- Dual safety straps that allow independent access to the Gastric and/or Jejunal ports
- External bolster minimizes contact and allows skin to breathe

**DEVICE MATERIALS**

The AMT Traditional G-JET<sup>®</sup> is composed of the following materials: Medical-grade silicone (96%) • Medical-grade thermoplastic (2%) • Medical-grade silicone pad print ink (1%) • Medical-grade, braided suture (1%) • Medical-grade thermoplastic (ENFit<sup>®</sup> Version Only) • Stainless steel for 16Fr and 18Fr devices

**DEVICE DURATION OF USE**

Device meets biocompatibility standard ISO 10993 and is approved for permanent contact (greater than 30 days).

**TRADITIONAL G-JET<sup>®</sup>**

Balloon feeding devices are meant to be periodically replaced for optimal performance, functionality, and cleanliness. Exact device longevity cannot be predicted. Device performance and functionality can degrade over time depending on usage and environmental conditions. Typical device longevity will vary for each patient depending on a number of factors, with typical device longevity ranging from 1-9 months. Some factors that can lead to reduced longevity include: gastric pH, diet of the patient, medications, balloon fill volume, trauma to the device, contact with sharp or abrasive objects, and overall tube care.

For optimal performance, it is recommended that the Traditional G-JET<sup>®</sup> device be changed at least every 3 months or as often as indicated by your healthcare professional. Proactive replacement of the device will help ensure optimal functionality and will help prevent unexpected device failure. If devices are failing or performance is degrading earlier than the typical range for device longevity, it is recommended that you speak with your healthcare professional regarding eliminating common factors that can lead to early device degradation. Also reference the **TROUBLESHOOTING** section for additional information regarding device performance issues.

**FEEDING AND DECOMPRESSION**

The following information explains how to properly feed into the small intestine and decompress the stomach.

**WARNING: THIS DEVICE HAS THE POTENTIAL TO MISCONNECT WITH SMALL BORE CONNECTORS OF OTHER HEALTHCARE APPLICATIONS. USE THIS DEVICE TO CONNECT ONLY TO COMPATIBLE ENTERAL DEVICES. DO NOT USE FOR NON-ENTERAL APPLICATIONS.**

**WARNING: WHEN USING A BOLUS STYLE CONNECTOR, THIS DEVICE HAS THE POTENTIAL TO MISCONNECT TO INTRAVENOUS, BREATHING SYSTEM, LIMB CUFF, AND TEMPERATURE CONNECTORS.**

**WARNING: TIGHTEN ONLY BY HAND. NEVER USE EXCESSIVE FORCE OR A TOOL TO TIGHTEN A ROTATING CONNECTOR. IMPROPER USE CAN LEAD TO CRACKING, LEAKAGE, OR OTHER FAILURE.**

**WARNING: ENSURE DEVICE IS CONNECTED TO AN ENTERAL PORT ONLY AND NOT TO AN IV SET.**

**WARNING: NEVER CONNECT THE JEJUNAL PORT TO SUCTION. DO NOT MEASURE RESIDUALS FROM THE JEJUNAL PORT.**

### JEJUNAL FEEDING

**CAUTION: Do not channel medications through the jejunal feeding port. This will clog the tube. To help avoid clogging, use a feeding pump to deliver formula through the lumen. Never allow formula to stand in the tube.**

1. Assemble the equipment: formula, feeding set, irrigating syringe, enteral feeding pump, and water for flushing the tube.
2. Wash your hands with soap and water. Rinse and dry them thoroughly.
3. Shake the formula and wipe the container tops before opening. If you do not use all the formula, cover the open container, record the date and time, and store it in the refrigerator. Throw away formula more than 24 hours old. Do not mix new formula with old formula. There is always a chance it will spoil.
4. Pour the formula into the feeding bag.
5. Inject the prescribed amount of water into the Traditional G-JET<sup>®</sup> by using the 35 ml irrigating syringe from your kit.
6. Connect the feeding set tubing to the enteral feeding pump. Follow the manufacturer's directions to set the pump.
7. If the feed set tubing has a clamp, open it.
8. Start the feeding pump.
9. Flush the jejunal port EVERY SIX HOURS by injecting the prescribed amount of water. **DO NOT USE FORCE.**

**NOTE:** If you are simultaneously decompressing the stomach while feeding and you see formula in the gastric drainage, stop the feeding and notify your specialist.

10. When the feeding is complete, clamp the feeding set tubing and disconnect the feeding set tubing from the feeding tube.

11. Flush the feeding tube thoroughly by injecting the prescribed amount of water.

**CAUTION: Do not use continuous or high intermittent suction. High pressure could collapse the tube or injure the stomach tissue and cause bleeding.**

### GASTRIC DECOMPRESSION

Some specialists recommend decompressing the stomach (letting the air or fluid out) before or after feeding. FOLLOW YOUR SPECIALIST'S INSTRUCTIONS. The gastric port may be left open for gravity drainage or connected to low intermittent suction. This allows stomach contents and/or gas to escape. Flush the gastric port EVERY SIX HOURS by injecting the prescribed amount of water. Do not use continuous or high intermittent suction. High pressure could collapse the tube or injure the stomach tissue and cause bleeding.

## TRADITIONAL G-JET<sup>®</sup> DEVICE – GENERAL USE AND CARE

### GENERAL FLUSHING GUIDELINES

Proper tube flushing is the best way to avoid clogging and maintain tube patency. The following are guidelines to avoid clogging and maintain tube patency.

- Flush the feeding tube with water every 4-6 hours during continuous feeding, anytime the feeding is interrupted, before and after every intermittent feeding, or at least every 8 hours if the tube is not being used.
- Use a 30 to 60 ml syringe. Do not use smaller size syringes as this can increase pressure on the tube and potentially rupture smaller tubes.
- Use room temperature tap water for tube flushing. Sterile water may be appropriate where the quality of municipal water supplies is of concern. The amount of water will depend on the patient's needs, clinical condition, and type of tube, but the average volume ranges from 10 to 50 ml for adults, and 3 to 10 ml for infants. Hydration status also influences the volume used for flushing feeding tubes. In many cases, increasing the flushing volume can avoid the need for supplemental intravenous fluid. However, individuals with renal failure and other fluid restrictions should receive the minimum flushing volume necessary to maintain patency.
- Flush the feeding tube before and after medication channeling and between medications. This will prevent the medication from interacting with formula and potentially causing the tube to clog.
- Use liquid medication when possible and consult the pharmacist to determine if it is safe to crush solid medication and to mix with water. If safe, pulverize the solid medication into a fine powder form and dissolve the powder in warm water before channeling through the feeding tube. Never crush enteric-coated medication or mix medication with formula.
- Avoid using acidic irrigants such as cranberry juice and cola beverages to flush feeding tubes as the acidic quality when combined with formula proteins may actually contribute to tube clogging.
- Do not use excessive force to flush the tube. Excessive force can perforate the tube and can cause injury to the gastrointestinal tract.
- Document the time and amount of water used in the patient's record. This will enable all caregivers to monitor the patient's needs more accurately.

### MEDICATION CHANNELING

Use liquid medication when possible and consult the pharmacist to determine if it is safe to crush solid medication and mix with water. If safe, pulverize the solid medication into a fine powder form and dissolve the powder in water before channeling through the feeding tube. Never crush enteric coated medication or mix medication with formula. Using a syringe, flush the tube with the prescribed amount of water.

### TUBING BLOCKAGE

Obstructions are usually caused by formula accumulation inside the tube. Body fluids and medications can also clog the tube. Prevent this by thoroughly flushing the tube EVERY SIX HOURS by injecting the prescribed amount of water. If the flow is sluggish, gently inject the water, and let it flow back out of the tube until the obstruction has cleared. When you feel too much resistance and cannot inject water at all, the tube may need to be repositioned or replaced because it may be kinked. High pressure can rupture the tube. Do not insert foreign objects into the tube. This may puncture the tube. If none of these measures work, call your specialist.

Tube occlusion is generally caused by: Poor flushing techniques • Failure to flush after measurement of gastric residuals • Inappropriate channeling of medication • Pill fragments • Viscous medications • Thick formulas, such as concentrated or enriched formulas that are generally thicker and more likely to obstruct tubes • Formula contamination that leads to coagulation • Reflux of gastric or intestinal contents up the tube

### TO UNCLOG A TUBE

**WARNING: DO NOT USE EXCESSIVE FORCE OR PRESSURE WHEN CHANNELING FEED OR MEDICATION, OR TO ATTEMPT TO CLEAR A CLOG IN THE TUBING. THIS CAN CAUSE THE TUBING TO RUPTURE OR LEAD TO DAMAGE OF THE TUBING SUPPORT STRUCTURE. IF CLOG CANNOT BE CLEARED OR TUBING IS CLOGGING ON A FREQUENT BASIS, THIS MAY INDICATE THAT THE DEVICE NEEDS TO BE REPLACED. FAILURE TO FOLLOW THIS WARNING MAY RESULT IN DEVICE FAILURE AND/OR SERIOUS PATIENT INJURY (E.G. BOWEL IRRITATION OR PERFORATION).**

1. Make sure that the feeding tube is not kinked or clamped off.
2. If the clog is visible above the skin surface, gently massage or milk the tube between fingers to break up the clog.
3. Next, place a syringe filled with warm water into the appropriate adaptor or lumen of the tube and gently pull back on, and then depress the plunger to dislodge the clog.
4. If the clog remains, repeat step #3. Gentle suction alternating with syringe pressure will relieve most obstructions.
5. If this fails, consult with the physician. Do not use cranberry juice, cola drinks, meat tenderizer or chymotrypsin, as they can actually cause clogs or create adverse reactions in some patients. If the clog is stubborn and cannot be removed, the tube will have to be replaced.

**GRANULATION TISSUE**

A small amount of extra tissue around the stoma site is simply the body's attempt to heal the wound. This is normal. It is also normal for the extra tissue to secrete a small amount of mucus. If the tissue bleeds excessively, drains, or interferes with the sliding external bolster, call your specialist.

**UNUSUAL SYMPTOMS**

Consult your health care practitioner if any of the following symptoms occur: Abdominal pain • Abdominal discomfort • Abdominal tenderness • Abdominal distension • Dizziness or fainting • Unexplained fever • Unusual amount of bleeding through or around the tube.

**DIARRHEA**

- The flow rate may be too fast. Decrease the rate, then increase it slowly until you reach the prescribed flow rate. Diluting the formula with water may help slow the diarrhea, but check with your specialist first.
- Prepare fresh formula for each feeding. Discard unused formula after 24 hours.
- Keep your equipment clean. Prepare only enough formula for eight hours. Wash, rinse and dry the equipment after each use.
- If the diarrhea continues for more than three days despite your efforts, call your specialist.
- Never attempt to bolus feed through the jejunal lumen.

**STOMA CARE**

Cleanse the site with normal saline three times daily to remove the small amount of mucus that normally accumulates around the stoma. If the mucus dries, it may irritate the skin. Applying a small amount of water may loosen the dried material. After the stoma heals, a thorough cleansing with soap and water is best. A dressing is not necessary, and contributes to skin maceration and infection. Unless recommended by your specialist, avoid medicated ointments or powders. DO NOT use oil or petroleum based products. DO NOT ROTATE THE TUBE. Jejunal tubes must NOT be rotated because they will kink. Should the tube become kinked, call your health care professional.

**BALLOON MAINTENANCE**

It is recommended that balloon volume is checked at least every two weeks, or as recommended by your healthcare professional.

**NOTE:** Refill the balloon using sterile or distilled water, not air or saline. Saline can crystallize and clog the balloon valve or lumen, and air may seep out and cause the balloon to collapse. Be sure to use the recommended amount of water as over-inflation can obstruct the lumen or decrease balloon life and under-inflation will not secure the tube properly.

**Healthcare Professionals:** Insert a Luer slip syringe into the balloon inflation port (**Fig. A-3**) and withdraw the fluid while holding the tube in place. Compare the amount of water in the syringe to the amount recommended or the amount initially prescribed and documented in the patient record. If the amount is less than recommended or prescribed, refill the balloon with the water initially removed, then draw up and add the amount needed to bring the balloon volume up to the recommended and prescribed amount of water. Be aware as you deflate the balloon there may be some gastric contents that can leak from around the tube. Document the fluid volume, the amount of volume to be replaced (if any), the date and time. Wait 10-20 minutes and repeat the procedure. The balloon is leaking if it has lost fluid, and the tube should be replaced. A deflated or ruptured balloon could cause the tube to dislodge or be displaced. If the balloon is ruptured, it will need to be replaced. Secure the tube into position using tape, then follow facility protocol and/or call the physician for instructions.

**End Users:** It is recommended that balloon volume is checked at least every two weeks, or as recommended by your healthcare professional. While holding the external bolster in place, insert a luer slip syringe into the balloon inflation port (**Fig. A-3**). Withdraw the water from the balloon. The amount of water you withdraw should equal the prescribed amount. If you withdraw less water than the prescribed amount, replace the amount recommended by your specialist. Never add water before you withdraw all of the water from the balloon. Note: The use of a luer-lock syringe may make this more difficult. A luer-slip syringe is recommended.

**DAILY CARE & MAINTENANCE CHECKLIST**

- **Assess the patient**
  - Assess the patient for any signs of pain, pressure or discomfort, warmth, rashes, purulent or gastrointestinal drainage.
  - Assess the patient for any signs of pressure necrosis, skin breakdown or hypergranulation tissue.
- **Clean the stoma site**
  - Use warm water and mild soap.
  - Use a circular motion moving from the tube outwards. Clean sutures, external bolsters and any stabilizing devices using a cotton-tipped applicator.
  - Rinse thoroughly and dry well.
- **Assess the tube**
  - Assess the tube for any abnormalities such as damage, clogging or abnormal discoloration.
- **Clean the feeding tube**
  - Use warm water and mild soap being careful not to pull or manipulate the tube excessively.
  - Rinse thoroughly, dry well.
- **Clean the jejunal, gastric and balloon ports**
  - Use a cotton tip applicator or soft cloth to remove all residual formula and medication.
- **Do not rotate the external bolster**
  - This will cause the tube to kink and possibly lose position.
- **Verify placement of the external bolster**
  - Verify that the sliding external bolster rests 2-3 mm above the skin.
- **Flush the feeding tube**
  - Flush the feeding tube with water every 4-6 hours during continuous feeding, anytime the feeding is interrupted, or at least every 8 hours if the tube is not being used.
  - Flush the feeding tube after checking gastric residuals. Flush the feeding tube before and after medication channeling. Avoid using acidic irrigants such as cranberry juice and cola beverages to flush feeding tubes.

**DEVICE REPLACEMENT**

**WARNING: REMOVAL IS RECOMMENDED TO BE PERFORMED BY A QUALIFIED CLINICIAN. SPONTANEOUS CLOSURE OF STOMA MAY OCCUR AS EARLY AS TWENTY-FOUR (24) HOURS AFTER REMOVAL. INSERT A NEW DEVICE IF ENTERAL FEEDING BY THIS ROUTE IS STILL INTENDED. IF CLOSURE IS DESIRED, APPLY A DRESSING OVER THE STOMA SITE.**

1. The Traditional G-JET® should be replaced periodically for optimal performance. Frequent performance checks are recommended. Clogging and/or reduced flow are indicators of diminished performance.
2. The Traditional G-JET® should only be used by or under the supervision of personnel trained in percutaneous gastrostomy tube placement. A thorough understanding of the technical principles, clinical applications, and risks associated with percutaneous gastrostomy tube placement is recommended prior to using this device.
3. Please contact a professional health care giver or physician for replacement of the Traditional G-JET®.

**NOTE:** The device can be disposed of by following local disposal guidelines, facility protocol or through conventional waste.

## MRI SAFETY INFORMATION

Non-clinical testing has demonstrated that the Traditional G-JET® feeding tube is MR Conditional. It can be scanned safely under the following conditions:

- Static magnetic field of 1.5-Tesla (1.5T) or 3.0-Tesla (3.0T).
- Spatial gradient field of up to:
  - 9,570 G/cm (95.70 T/m) for 1.5T systems
  - 5,720 G/cm (57.20 T/m) for 3.0T systems
- Maximum whole body averaged specific absorption rate (SAR) of:
  - 4.0 W/kg for 15 minutes of scanning in Normal Operating Mode at 1.5T.
  - 4.0 W/kg for 15 minutes of scanning in Normal Operating Mode at 3.0T.

**1.5 RF heating:** In non-clinical testing with body coil excitation, the Traditional G-JET® feeding tube produced a temperature rise of less than 1.0°C at a maximum whole body averaged specific absorption rate (SAR) of 4.0 W/kg, as assessed by calorimetry for 15 minutes of scanning in a 1.5T Siemens Espree (MRC30732) MR scanner with SYNGO MR B17 Software.

**3.0 RF heating:** In non-clinical testing with body coil excitation, the Traditional G-JET® feeding tube produced a temperature rise of less than 1.0°C at a maximum whole body averaged specific absorption rate (SAR) of 4.0 W/kg, as assessed by calorimetry for 15 minutes of scanning in a 3.0T Siemens Trio (MRC20587) MR scanner with SYNGO MR A30 4VA30A software.

**3.0T MR Artifact:** In testing using a 3.0T system with gradient-echo sequencing, the shape of the image artifact follows the approximate contour of the device and extends radially up to the 2.0 cm from the implant.

**NOTE:** A magnetic resonance environment safety card has been provided in this kit. Please keep this card with you for your records. This card can be given to your radiologist when undergoing an MRI scan.

## TROUBLESHOOTING

Long-term device performance and functionality are dependent on proper device usage per the instructions and varying usage and environmental factors. While it is expected that you will use your feeding device without any issues, unexpected device problems can sometimes occur. The following section covers a number of performance or functionality related items and how to help prevent these types of occurrences.

**Gastric leakage:** Gastric leakage may occur if the sliding external bolster is not properly adjusted or if the balloon is not against the stomach wall. Adjust the external bolster so that it sits 1-2 mm from the skin. Once the bolster is adjusted, check the water volume in the balloon. Do not add more water to the balloon before first withdrawing all of the water from the balloon. Fill the balloon according to the balloon inflation range in Table 1.

**Leaking of balloon volume:** If balloon is deflating, fully deflate the balloon. Compare the amount of water in the syringe to the amount recommended or the amount initially prescribed and documented in the patient record. Inflate balloon with recommended fill-volume. Do not exceed maximum fill volume. Wait 10-20 minutes and repeat the procedure. The balloon is leaking if it has lost fluid, and the tube should be replaced. Only access the balloon inflation port for balloon inflation and deflation. Other uses can result in balloon leakage or device failure. **NOTE:** The balloon is made of a semi-permeable material and can lose a small amount of fill-volume over time depending on environment and usage conditions.

**A tear has formed:** Tears can occur due to contact with a sharp or abrasive object, excessive force, or excessive pressure. Due to the soft, comfortable nature of the material the device is made from, small tears can quickly lead to large tears or device failure. If a tear is noticed on the device, consider replacing device and check for any sources of tension, force, or sharpness that may be leading to the tears occurring.

**Balloon inflation valve leakage:** Leakage from this valve typically occurs due to residue in the fill-valve. Always use a clean syringe when inflating the balloon and only inflate water through the valve. Leakage can also occur due to valve becoming stuck if syringe is inserted too hard into valve. Insert the syringe into the balloon inflation port if leakage is occurring to attempt to reset valve. Several attempts may be necessary before valve resets.

**Foul smell coming from the device:** Foul smells can occur due to not properly flushing the device after each use, infection, or other growth forming inside of the device. If a foul smell is noted coming from the device, device should be flushed and stoma site should be gently cleaned with soap and warm water. If foul smell does not go away, it is recommended that you contact your healthcare professional.

**Plug will not stay closed:** Ensure that the plug is being firmly and fully pressed or rotated without excessive force. If the plug is not staying closed, check the plug and feed-port area for any excess residue build-up. Clean excess residue build-up with cloth and warm water.

**Balloon will not inflate or deflate:** Inflation and deflation problems typically occur due to residue blocking the fill-lumen. Always use a clean syringe when inflating the balloon. In some rare occurrences, fungus can grow inside of the tubing and block the fill-lumen. Fungus growth can occur depending on patient environment and feeds/medications being channeled through the device. If balloon will not deflate, contact your healthcare professional for removal. If deflation problem is the result of fungal growth, eliminating source of fungal growth or anti-fungal medication may be required.

**Balloon is misshaped:** Be sure to inflate and inspect balloon prior to placement to check balloon symmetry. Balloons becoming misshaped typically occur due to excessive force or tension on the device (pulling device out of stoma while balloon is inflated). A device that is placed too close to the pylorus can result in a misshaped balloon and blockage of the pylorus. Balloons can be gently massaged with fingers back into symmetry if not excessively misshaped. Device should be replaced if balloon is found excessively misshaped.

**Tubing has reduced flow or has become clogged:** Tubing can become blocked due to not properly flushing after each use, use of thick or improperly crushed medications, use of thick feeds/formulas, gastric reflux, and/or fungus growth. If clogged, refer to the **GENERAL USE AND CARE** section for instructions on how to unclog the device. If clog cannot be removed, device may need to be replaced. **WARNING: DO NOT USE EXCESSIVE FORCE OR PRESSURE WHEN CHANNELING FEED OR MEDICATION, OR TO ATTEMPT TO CLEAR A CLOG IN THE TUBING. THIS CAN CAUSE THE TUBING TO RUPTURE OR LEAD TO DAMAGE OF THE TUBING SUPPORT STRUCTURE. IF CLOG CANNOT BE CLEARED OR TUBING IS CLOGGING ON A FREQUENT BASIS, THIS MAY INDICATE THAT THE DEVICE NEEDS TO BE REPLACED. FAILURE TO FOLLOW THIS WARNING MAY RESULT IN DEVICE FAILURE AND/OR SERIOUS PATIENT INJURY (E.G. BOWEL IRRITATION OR PERFORATION).**

**Balloon failure:** Early balloon failure can occur due to a number of patient or environmental factors, including but not limited to: gastric pH, diet, certain medications, improper balloon fill volume, placement of device, improperly channeling food/medication into balloon port, trauma, contact with a sharp or abrasive material, and overall care of the device.

**Device has become discolored:** The device can become discolored over days to months of use. This is normal depending on the types of feeds and medications being used with the device.

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ENFit® is a registered trademark of Global Enteral Device Supplier Association, Inc.
**Applied Medical Technology, Inc.**