# ENFit<sup>®</sup> Transition Adapters

## Instructions for Use

CE  
2862

- FR** Adaptateurs de transition ENFit<sup>®</sup>. Mode d'emploi.
- ES** Adaptadores de transición ENFit<sup>®</sup>. Instrucciones de uso.
- DE** ENFit<sup>®</sup> -Übergangsadapter. Gebrauchsanweisung.
- PT-BR** Adaptadores de Transição ENFit<sup>®</sup>. Instruções de utilização.
- IT** Adattatore da ENFit<sup>®</sup>. Istruzioni per l'uso.
- NO** ENFit<sup>®</sup> overgangsadaptere. Brukerveiledning.
- SV** ENFit<sup>®</sup> Övergångsadaptrar. Bruksanvisning.
- FI** ENFit<sup>®</sup>-siirtoadapterit. Käyttöohjeet.
- NL** ENFit<sup>®</sup> Transitie-adapters. Instructies voor gebruik.
- DA** ENFit<sup>®</sup>-overgangsadaptere. Brugsvejledning.
- AR** تعليمات الاستخدام ENFit<sup>®</sup>, المهايئات الانتقالية
- BG** Преходни адаптери ENFit<sup>®</sup>. Инструкции за употреба.
- CS** Přechodové adaptéry ENFit<sup>®</sup>. Návod k použití.
- EL** ENFit<sup>®</sup> Προσαρμογέας μετάδοσης. Οδηγίες χρήσης.
- ET** Üleminekuadapterid ENFit<sup>®</sup>. Kasutusjuhised.
- HU** ENFit<sup>®</sup> átmeneti adapterek. Használati utasítások.
- JA** ENFit<sup>®</sup> トランジション・アダプタ. 使用説明書.
- PL** Przejściówki ENFit<sup>®</sup>. Instrukcja stosowania.
- RO** Adaptoare de tranziție ENFit<sup>®</sup>. Instrucțiuni de utilizare.
- SL** Prehodni prilagodilniki ENFit<sup>®</sup>. Navodila za uporabo.
- TR** ENFit<sup>®</sup> Geçiş Adaptörleri. Kullanma Talimatları.
- LT** „ENFit“ perjungimo adapteriai. Naudojimo instrukcija.
- PT** Adaptadores de Transição ENFit<sup>®</sup>. Instruções para a Utilização.
- SK** Prechodné adaptéry ENFit<sup>®</sup>. Návod na používanie.
- LV** ENFit<sup>®</sup> pārejas adapter. Lietošanas instrukcijas.
- MS** Penyesuai Peralihan ENFit<sup>®</sup>. Arahan penggunaan.

R<sub>x</sub> ONLY

C4510-E 03/2024

Made in the U.S.A.

# ENFit® Transition Adapters

## Instructions for Use

### INTENDED USE

ENFit® Transition Adapters are intended to be used as effective adapters that facilitate connections between legacy and ENFit® enteral connectors. ENFit® Transition Adapters are intended to be used by clinicians and caregivers for infant, child, adolescent, adult, and elderly patients. The ENFit® Transition Adapters aid patients with a wide array of medical conditions that make maintaining an adequate nutritional state a challenge.

### INDICATIONS FOR USE

The AMT Enteral Transition Adapters are intended to facilitate enteral specific connections between AAMI/CN3(PS) compliant connectors and non ISO 80369-1 compliant legacy enteral connectors.

### CONTRAINDICATIONS

At this time there are no contraindications for use of the ENFit® Transition Adapters.

### COMPLICATIONS

At this time there are no complications associated with use of the ENFit® Transition Adapters.

**NOTE:** Please contact AMT, our Authorized Representative (EC Rep), and/or the competent authority of the member state in which you are established if a serious incident has occurred in relation to the device.

### CLINICAL BENEFITS

Clinical benefits to be expected when using the ENFit® Transition Adapters include but are not limited to: Provides a connection method between feeding sets, extension sets and syringes with ENFit® connectors and enteral connections with non-ENFit® compliant connecting ends

### PERFORMANCE CHARACTERISTICS

Performance characteristics of the ENFit® Transition Adapters include but are not limited to: Provide connections between giving sets and syringes with ISO 80369-3 ENFit® connectors and enteral connections with non-ISO 80369-3 compliant connecting ends • Wide variety of ENFit® Transition Adapters • Forward and backwards compatibility - the ability to connect current market to ENFit® and connect ENFit® to current market • Glow Green Technology provides better visibility in the dark for easier nighttime feeds • Easily removed for cleaning between uses

### INSTRUCTIONS FOR USE

**WARNING: THIS PACKAGE CONTAINS SMALL PARTS THAT ARE A POTENTIAL CHOKING HAZARD IF NOT USED AS INTENDED. USE UNDER APPROPRIATE ADULT SUPERVISION.**

**WARNING: THIS TRANSITION ADAPTER DOES NOT PROTECT AGAINST INADVERTENT MISCONNECTIONS WITH CONNECTORS OF NON-ENTERAL HEALTHCARE DEVICES. FOR ENTERAL USE ONLY.**

**WARNING: IF USING A LEGACY-STYLE (NON-ENFIT®) CONNECTOR, THIS DEVICE HAS THE POTENTIAL TO MISCONNECT TO THE FOLLOWING SYSTEMS: BREATHING, CONES AND SOCKETS OF ANESTHETIC AND RESPIRATORY EQUIPMENT, INTRAVENOUS, LIMB CUFF, NEURAXIAL CONNECTORS, NIPPLES OF RESPIRATORY THERAPY EQUIPMENT, URINARY, AND TEMPERATURE SENSOR CONNECTORS OF RESPIRATORY HUMIDIFYING EQUIPMENT.**

1. Inspect contents for damage. If damaged, do not use. Obtain another package.
2. Use the included adapter to connect a AAMI/CN3(PS) compliant connector to non ISO 80369-1 compliant legacy enteral connector. Screw the adapter clockwise to lock the ENFit® connectors together. Catheter and luer connectors are compression fit.

**WARNING: IF THE ADAPTER IS NOT PROPERLY PLACED AND LOCKED, LEAKAGE MAY OCCUR. DURING CONNECTION, APPLY ENOUGH FORCE TO THE ADAPTERS TO PROVIDE A SECURE CONNECTION. DO NOT OVER-TIGHTEN THE ADAPTER TO PREVENT DAMAGE TO THE CONNECTORS.**

3. The enteral adapter should be removed after each use. Unscrew the ENFit® adapter with a counter-clockwise turn. Catheter and luer connectors are compression fit and can be removed using moderate tension.
4. Clean the adapter with warm soapy water and thoroughly rinse and dry after each use. Store in a clean, dry place. Never use a dishwasher for cleaning.

#### DEVICE LONGEVITY


Always remove the transition adapter from the legacy adapter after each use. The transition adapter is re-usable. Depending on usage (frequency, environment, etc.), periodic replacement of the adapter may be necessary for optimal functionality. Replace adapter if any leakage or cracks occur.

**NOTE:** The device can be disposed of by following local disposal guidelines, facility protocol or through conventional waste.

**NOTE:** Refer to your feeding device Directions For Use for the circumstances in which you should consult a healthcare professional.

#### THANK YOU!

Thank you for choosing AMT. For additional help and information regarding the use of our device, feel free to contact AMT with the contact information on the back of the instructions for use. We are happy to hear your thoughts and help with your concerns and questions.

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© 2024 Applied Medical Technology, Inc.
**Applied Medical Technology, Inc.**