Introducer/Stiffener for Initial Placement

CE
2862

- ES** Introducitor-guía para la colocación inicial
DE Inserter/Versteifer für die anfängliche Platzierung
IT Cannula di introduzione/elemento di irrigidimento per il posizionamento iniziale
FR Introducteur/Raidisseur pour la mise en place initiale
NO Introduser/oppstrammer for innledende plassering
SV Introducer/Ledare för initial placering
NL Introducer/versteviger voor initiële plaatsing
DA Introducer/afstiver til førstegangsanlæggelse
FI Sisäänviejä/tuki ensiasennusta varten
AR منخال / دعامة الموضعة الابتدائية
JA 初期配置用インテロデューサ/スティフナ
PL Introduktor/usztywniacz do zakładania początkowego
BG Устройство за въвеждане/Арматура за първоначално поставяне
CS Zaváděč/výztuha k úvodnímu umístění
ET Sisesti/jäigastaja esmapaigalduseks
EL Εισαγωγικός εισαγωγέας για αρχική τοποθέτηση
HU Bevezető/merevítő kezdeti behelyezéshez
LV Sākotnējās ievietošanas ievadītājs/cietinātājs
PT-BR Introduutor /Enrijecedor para colocação inicial
SL Uvajalo/ojačevalo za začetno vstavitv
TR İlk Yerleştirme İçin İntrodüser/Sertleştirici
RO Dispozitiv de introducere/fixare pentru plasarea inițială
LT Introduiseris / standiklis pirminiam įdėjimui
PT Introduutor/Enrijecedeor para Colocação Inicial
SK Zavádzač/výstuž pre počiatočné zavedenie

R_X ONLY

DEHP

C4369-E 07/2025

Made in U.S.A.

STERILE EO

DIRECTIONS FOR USE

CAUTION: Federal (USA) law restricts this device to sale, distribution, and use by or on the order of a physician.

DESCRIPTION:

Optional Introducer/Stiffener makes insertion of device easier and is guidewire compatible.

INTENDED USE:

The Introducer/Stiffener for Initial Placement is an effective tool used to stiffen the MINI™ Button for easier placement. The device is especially useful in patients with stoma tracts that are long, tight or angled. The Introducer/Stiffener is intended to be used by clinicians for patients of all ages (infant to elderly).

INDICATIONS FOR USE:

The Introducer is indicated to assist in the placement of AMT MINI™ Button feeding devices by providing temporary feeding tube stability during insertion through well-established patient stoma tract.

CONTRAINDICATIONS FOR USE:

1. Lack of adherence of the stomach to the abdominal wall
2. Lack of a well-established gastrostomy site
3. Evidence of infection
4. Uncertainty as to gastrostomy tract direction and length (abdominal wall thickness)

COMPLICATIONS:

At this time there are no complications associated with the use of the Introducer/Stiffener for Initial Placement.

NOTE: Please contact AMT, our Authorized Representative (EC Rep), and/or the competent authority of the member state in which you are established if a serious incident has occurred in relation to the device.

CLINICAL BENEFITS:

Clinical Benefits to be expected when using the Introducer/Stiffener for Initial Placement include but are not limited to:

- Provides stiffening support to the MINI™ Button for easier placement; especially true for patients with long, tight, or angled stoma tracts
- Increased ease of placement reduces pain from replacing the device while ensuring the patient can receive proper nutrition and medication
- Guidewire compatibility helps to ensure the MINI™ Button is in the stomach and not the abdominal wall

PERFORMANCE CHARACTERISTICS:

Performance Characteristics of the Introducer/Stiffener for Initial Placement include but are not limited to:

- Provides ease of placement of MINI™ Button
- Compatible with devices 12F and larger (4-6000, 4-6200)
- Compatible with devices 10F and larger (4-6010)
- Compatible with devices up to 6.5 cm in length
- Guidewire compatible

1. DEVICE PREPARATION

CAUTION: Prior to placement, please inspect all contents of the kit for damage. If the package is damaged or sterile barrier is breached, do not use the product.

- During device preparation, lubricate the tip portion of the enclosed Introducer/Stiffener using water-soluble lubricant.
- Insert the Introducer/Stiffener into the device opening until the tip is flush with the tip of the device and necessary stiffening is achieved.
- Lubricate the outer surface of the feeding device using water-soluble lubricant.

2. DEVICE PLACEMENT

- The included Introducer/Stiffener is guidewire compatible. During placement of a feeding device with a guidewire already in place, pass the Introducer/Stiffener over the guidewire.
- Follow instructions for feeding device for remainder of device placement.

3. STIFFENER REMOVAL

- After feeding device is placed, slowly withdraw the Introducer/Stiffener and guidewire from the device.

CAUTION: Care must be taken when withdrawing the Introducer/Stiffener to ensure that the entire device is not also removed.

WARNING: THIS DEVICE IS INTENDED FOR SINGLE USE. DO NOT REUSE, REPROCESS OR RE-STERILIZE THIS MEDICAL DEVICE. DOING SO MAY COMPROMISE BIOCOMPATIBILITY CHARACTERISTICS, DEVICE PERFORMANCE AND/OR MATERIAL INTEGRITY; ANY OF WHICH MAY RESULT IN POTENTIAL PATIENT INJURY, ILLNESS AND/OR DEATH.

NOTE: The device can be disposed of by following local disposal guidelines, facility protocol or through conventional waste.

THANK YOU!

Thank you for choosing AMT. For additional help and information regarding the use of our device, feel free to contact AMT with the contact information on the back of the instructions for use. We are happy to hear your thoughts and help with your concerns and questions.

 Rx ONLY Prescription Only	 Sterilized Using Ethylene Oxide
 Not made with natural rubber latex	 Not made with DEHP (Di(2-ethylhexyl) Phthalate)
 Do not resterilize	 Do not use if package is damaged and consult instructions for use
 Single Use Only	 Medical Device

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Applied Medical Technology, Inc.