Patient Information Leaflet

MiniACE® Button

Low Profile Antegrade Continence Enema Button

ENGLISH

What is in this leaflet?

This leaflet will answer some questions about the MiniACE® Button. AMT has provided this information as an educational resource tool. This is not intended as a substitute for professional medical care. Your FIRST source of information should be your healthcare provider.

What is the MiniACE® Button?

The MiniACE® Button is a type of antegrade continence enema tube that is inserted directly into the colon through an opening known as a stoma. The device includes a balloon that is filled with distilled or sterile water to help secure it inside the colon and prevent it from falling out, along with an external bolster to help maintain tube position. The device is available in various French sizes and lengths to accommodate the needs of patients.

The MiniACE® Button is made from medical-grade silicone (81%), medical-grade thermoplastic (16%), stainless steel spring (2%), and medical-grade silicone pad print ink (1%). There are no manufacturing residuals that could pose a risk to the patient.

What is the MiniACE® Button Used for?

The MiniACE® Button is intended to be used as an effective device for channeling irrigation fluid through a secured (initial placement) or formed (replacement) stoma. The MiniACE® Button is intended to instill fluids through a stoma into the colon to promote evacuation of the contents of the lower bowel through the anus and is intended to be an aid in the management of fecal incontinence. The catheter is placed and maintained in a percutaneously prepared opening, such as a cecostomy or appendicostomy. The device is intended to be used in children and adults.

How do you use the MiniACE® Button after it has been placed?

Irrigation: Ensure the irrigation set clamp is closed and attach the irrigation set to the button by lining up the dark line on the irrigation set connector with the dark line of the interlock on the button. Fully press the irrigation set connector into the button. Turn ¾ to the right (clock-wise) to lock the irrigation set in place. Attach the opposite end of the irrigation set to the connector being used. Never use excessive force or a tool to tighten a rotating connector. Improper use can lead to cracking, leakage, or other failure. Once connected, open the clamp to allow flow.


Figure 2: Attaching the irrigation set

When irrigation is complete, flush with 5-10 ml water (if not using water for irrigation). Close the clamp to prevent leakage while removing the irrigation set. Remove the irrigation set by holding the Button in place, grasp the irrigation connector, twist ¾ turn to the left (counter clock-wise) to unlock, and gently remove the irrigation set from the Button. Snap the MiniACE® Button plug in place to keep the lumen clean.


Figure 1: MiniACE® Button

How do you take care of your MiniACE® Button after it has been placed?

Device Care: Check the device daily for any signs of damage or clogging. Clogging and/or reduced flow are indications of diminished performance. The stoma area should be kept clean and dry at all times. It is important to clean the stoma site every day. A cotton swab or terry cloth can be used to clean the skin around the device with mild soap and water. Clean the stoma as directed by your healthcare professional. The device should be rotated daily for site hygiene. Always check the stoma site for redness, pain/soreness, swelling, or any drainage. In case of fever, abdominal distention, infection, blockage, or tissue necrosis, patients should see their physician immediately. If any of these symptoms are observed, contact your healthcare professional for advice.

Clogging: First check to make sure that the tube is not kinked or clamped anywhere. If there is a visible clog in the tubing, attempt to massage the device to break up the clog. Connect a syringe to an irrigation set and attach into the interlock connector. Fill the syringe with warm water and gently push and pull the syringe plunger to free the clog. It may take several cycles of pushing / pulling the plunger to clear the clog. If clog cannot be removed, contact your healthcare professional, as the tube may need to be replaced. Do not use excessive force or pressure to attempt to clear the clog. This can cause the tubing to rupture.

Balloon Care: It is recommended that the balloon volume is checked as instructed by your healthcare professional. Remove the water with a syringe and compare the amount removed to the recommended amount. Refill the balloon and, if needed, add additional water to meet the recommended amount. Wait 10-20 minutes and repeat. The balloon is leaking if it has lost fluid, and the tube should be replaced. If the balloon is damaged, secure the tube in place using tape, then call your healthcare professional for instructions.

Refill the balloon using sterile or distilled water, not air or saline. Saline can crystallize and clog the balloon valve or lumen, and air may seep out and cause the balloon to collapse. Be sure to use the recommended amount of water as over-inflation can obstruct the lumen or decrease balloon life and under-inflation will not secure the tube properly.

How long does the MiniACE® Button last?

Low-profile balloon cecostomy devices are meant to be periodically replaced for optimal performance and functionality. Device performance and functionality can degrade over time depending on usage and environmental conditions. Typical device longevity will vary for each patient depending on a number of factors, with typical device longevity ranging from 1-6 months. Some factors that can lead to reduced longevity include: intestinal pH, diet of the patient, medications, balloon fill volume, trauma to the device, contact with sharp or abrasive objects, incorrect stoma length measurement, and overall tube care.

For optimal performance, it is recommended that MiniACE® Button be changed at least every 3 months or as often as indicated by your healthcare professional. Proactive replacement of the device will help ensure optimal functionality and will help prevent unexpected failure.

What are the possible side effects of having the MiniACE® Button?

Potential complications when using the MiniACE® Button include but are not limited to: Peristomal pain • Abscess, wound infection and skin breakdown • Pressure necrosis • Hypergranulation tissue • Intraperitoneal leakage • Buried bumper syndrome • Peristomal leakage • Balloon failure or dislodgement • Tube clog • Colonic bleeding and/or ulcerations • Bowel volvulus

Please consult your healthcare professional if you experience any of the following: Fever, vomiting, or diarrhea • Skin around the stoma site is red, discolored or raw • Drainage around the stoma site is white, yellow, or green; drainage has an unpleasant odor • Crusting is noted at the stoma site • Large amount of tissue build-up (such as granulation tissue) • Swollen skin or tissue at the stoma site • Repetitive leakage of intestinal contents • Pain, bleeding, puss, or inflammation at stoma site • Device is no longer fitting properly • Device falls out and cannot be easily replaced • Abdominal distention

There is a risk that if the internal balloon fails, the tube could fall out. There could be leakage of intestinal content around the tube or leakage from the device. The tube could become clogged or have reduced flow. The interlock (irrigation port) may crack, leak, or separate from the device. The device may become discolored over days to months of use.

Is the MiniACE® Button MRI Compatible?

Non-clinical testing demonstrated that the MiniACE® Button is MR Conditional. A patient with this device can be scanned safely in an MR system under the following conditions:

- Static magnetic field of 1.5-Tesla and 3-Tesla, only
- Maximum spatial gradient magnetic field of 1,000-gauss/cm (10-T/m)
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning (i.e., per pulse sequence) in the Normal Operating Mode

Under the scan conditions defined, the MiniACE® Button is expected to produce a maximum temperature rise of 1.6°C after 15 minutes of continuous scanning (i.e., per pulse sequence).

In non-clinical testing, the image artifact caused by the MiniACE® Button extends approximately 20-mm from this device when imaged using a gradient echo pulse sequence and a 3-Tesla MR system.

How do you report adverse events?

- For Australian Users: Notice that any serious incident that occurs in relation to the device should be reported to Applied Medical Technology, Inc. and The Therapeutic Goods Administration (TGA) at <http://www.tga.gov.au>.
- Please contact AMT, our Authorized Representative (EC Rep), and/or the competent authority of the member state in which you are established if a serious incident has occurred in relation to the device.

Product Codes:

Length (cm)	MiniACE® Button Kits (Legacy)			MiniACE® Button Kits (ENFit®)		
0.8	ACE-1008	ACE-1208	ACE-1408	ACE-1008-I	ACE-1208-I	ACE-1408-I
1.0	ACE-1010	ACE-1210	ACE-1410	ACE-1010-I	ACE-1210-I	ACE-1410-I
1.2	ACE-1012	ACE-1212	ACE-1412	ACE-1012-I	ACE-1212-I	ACE-1412-I
1.5	ACE-1015	ACE-1215	ACE-1415	ACE-1015-I	ACE-1215-I	ACE-1415-I
1.7	ACE-1017	ACE-1217	ACE-1417	ACE-1017-I	ACE-1217-I	ACE-1417-I
2.0	ACE-1020	ACE-1220	ACE-1420	ACE-1020-I	ACE-1220-I	ACE-1420-I
2.3	ACE-1023	ACE-1223	ACE-1423	ACE-1023-I	ACE-1223-I	ACE-1423-I
2.5	ACE-1025	ACE-1225	ACE-1425	ACE-1025-I	ACE-1225-I	ACE-1425-I
2.7	ACE-1027	ACE-1227	ACE-1427	ACE-1027-I	ACE-1227-I	ACE-1427-I
3.0	ACE-1030	ACE-1230	ACE-1430	ACE-1030-I	ACE-1230-I	ACE-1430-I
3.5	ACE-1035	ACE-1235	ACE-1435	ACE-1035-I	ACE-1235-I	ACE-1435-I
4.0	ACE-1040	ACE-1240	ACE-1440	ACE-1040-I	ACE-1240-I	ACE-1440-I
4.4	ACE-1044	ACE-1244	ACE-1444	ACE-1044-I	ACE-1244-I	ACE-1444-I
5.0	ACE-1050	ACE-1250	ACE-1450	ACE-1050-I	ACE-1250-I	ACE-1450-I
5.5	ACE-1055	ACE-1255	ACE-1455	ACE-1055-I	ACE-1255-I	ACE-1455-I
6.0	ACE-1060	ACE-1260	ACE-1460	ACE-1060-I	ACE-1260-I	ACE-1460-I
6.5	ACE-1065	ACE-1265	ACE-1465	ACE-1065-I	ACE-1265-I	ACE-1465-I
7.0	ACE-1070	ACE-1270	ACE-1470	ACE-1070-I	ACE-1270-I	ACE-1470-I
8.0	ACE-1080	ACE-1280	ACE-1480	ACE-1080-I	ACE-1280-I	ACE-1480-I
9.0	ACE-1090	ACE-1290	ACE-1490	ACE-1090-I	ACE-1290-I	ACE-1490-I
10.0	ACE-10100	ACE-12100	ACE-14100	ACE-10100-I	ACE-12100-I	ACE-14100-I

Applied Medical Technology, Inc.

8006 Katherine Boulevard

Brecksville, OH 44141 USA

Toll Free: +1-800-869-7382 | Telephone: +1-440-717-4000

Website: www.appliedmedical.net | E-mail: cs@appliedmedical.net

© 2025 Applied Medical Technology, Inc. Rev C4731-D 07/2025

Patents: <http://www.AppliedMedical.net/Patents/> | Trademarks: <http://www.AppliedMedical.net/Trademarks/>

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