Patient Information Leaflet

MiniONE® Non-Balloon Button

Low Profile Gastrostomy Feeding Tube

ENGLISH

What is in this leaflet?

This leaflet will answer some questions about the MiniONE® Non-Balloon Button. AMT has provided this information as an educational resource tool. This is not intended as a substitute for professional medical care. Your FIRST source of information should be your healthcare provider.

What is the MiniONE® Non-Balloon Button?

The MiniONE® Non-Balloon Button is a type of gastrostomy feeding tube (G-tube) that is inserted directly into the stomach through an opening known as a stoma. The device includes an internal silicone retention bolster to help secure it inside the stomach and prevent it from falling out, along with an external bolster to help maintain tube position. The device is available in various French sizes and lengths to accommodate the needs of patients of all age groups.

The MiniONE® Non-Balloon Button is made from medical-grade silicone (92%), medical-grade thermoplastic (7%), and medical-grade silicone pad print ink (1%). There are no manufacturing residuals that could pose a risk to the patient.

What is a MiniONE® Non-Balloon Button Used for?

The MiniONE® Non-Balloon Button provides nutrition, medication, and decompression access into the stomach through a secured (initial placement) or formed (replacement) stoma. The device is indicated for use in patients who require long term feeding, are unable to tolerate oral feeding, who are at low risk for aspiration, require gastric decompression and/or medication delivery directly into the stomach. The device is intended for all age groups.

How do you use a MiniONE® Non-Balloon Button after it has been placed?

Feeding: Ensure the feeding set clamp is closed and attach the feeding set to the button by lining up the dark line on the feeding set connector with the dark line on the button. Fully press the feeding set connector into the button. Turn $\frac{3}{4}$ to the right (clock-wise) to lock the feeding set in place. Attach the opposite end of the feeding set to the connector being used. Never use excessive force or a tool to tighten a rotating connector. Improper use can lead to cracking, leakage, or other failure. Once connected, open the clamp to allow flow.

Figure 2: Attaching the feeding set

Medication Channeling: Liquid medication is preferred when available. If solid medication is required, consult with your physician regarding whether it is safe to crush the medication. If safe, medication should be crushed as fine as possible (into powder form) and dissolved into water before channeling medication through device. Channeling solid medication through the tubing that is not properly crushed can lead to blockage in the tubing. Never crush enteric coated medication or mix medication with formula. Flush device with water after channeling medication.


Figure 1: MiniONE® Non-Balloon Button

Decompression: Decompression can be performed if there is a need to check for residuals, or if the patient tends to retain air during feeding causing bloating and discomfort. To decompress, connect the feeding set to the interlock connector and drain stomach contents into a container. After decompression is complete, flush the feeding set with water and remove the feeding set from the device.

Flushing: Always use a MiniONE® feeding set to flush the MiniONE® Non-Balloon Button. Use room temperature water for tube flushing. The amount of water will depend on the patient's needs, clinical condition, and type of tube, but the average volume ranges from 10 to 50 ml for adults, and 3 to 10 ml for infants. Flush the feeding tube with water every 4-6 hours during continuous feeding, anytime the feeding is interrupted, before and after every intermittent feeding, or at least every 8 hours if the tube is not being used. Flush the feeding tube before and after medication channeling and between medications. Do not use excessive force to flush the tube. Excessive force can perforate the tube and can cause injury to the gastrointestinal tract.

How do you take care of your MiniONE® Non-Balloon Button after it has been placed?

Device Care: Check the device daily for any signs of damage or clogging. Clogging and/or reduced flow are indications of diminished performance. The stoma area should be kept clean and dry at all times. It is important to clean the stoma site every day. A cotton swab or terry cloth can be used to clean the skin around the device with mild soap and water. Clean the stoma as directed by your healthcare professional. The device should be rotated daily for site hygiene. Always check the stoma site for redness, pain/soreness, swelling, or any drainage. If any of these symptoms are observed, contact your healthcare professional for advice.

Clogging: First check to make sure that the tube is not kinked or clamped anywhere. If there is a visible clog in the tubing, attempt to massage the device to break up the clog. Connect a syringe to an extension set and attach into the interlock connector. Fill the syringe with warm water and gently push and pull the syringe plunger to free the clog. It may take several cycles of pushing / pulling the plunger to clear the clog. If clog cannot be removed, contact your healthcare professional, as the tube may need to be replaced. Do not use excessive force or pressure to attempt to clear the clog. This can cause the tubing to rupture.

How long does a MiniONE® Non-Balloon Button last?

Low-profile feeding devices are meant to be periodically replaced for optimal performance, functionality, and cleanliness. Device performance and functionality can degrade over time depending on usage and environmental conditions. Typical device longevity will vary for each patient depending on a number of factors, with typical device longevity ranging from 3-9 months. Some factors that can lead to reduced longevity include: gastric pH, diet of the patient, medications, trauma to the device, contact with sharp or abrasive objects, incorrect stoma length measurement, and overall tube care.

For optimal performance, it is recommended that MiniONE® Non-Balloon Button device be changed every 6 months or as often as indicated by your healthcare professional. Proactive replacement of the device will help ensure optimal functionality and will help prevent unexpected device failure.

What are the possible side effects of having a MiniONE® Non-Balloon Button?

Potential complications when using the MiniONE® Non-Balloon Button include but are not limited to: Nausea, vomiting, abdominal bloating or diarrhea • Aspiration • Peristomal pain • Abscess, wound infection and skin breakdown • Pressure necrosis • Hypergranulation tissue • Intraperitoneal leakage • Buried bumper syndrome • Peristomal leakage • Tube clog • Gastrointestinal bleeding and/or ulcerations • Ileus or gastroparesis • Bowel and gastric volvulus • Peritonitis • Gastrocolic fistula • Sepsis

Please consult your healthcare professional if you experience any of the following: Fever, vomiting, or diarrhea • Skin around the stoma site is red, discolored or raw • Drainage around the stoma site is white, yellow, or green; drainage has an unpleasant odor • Crusting is noted at the stoma site • Large amount of tissue build-up (such as granulation tissue) • Swollen skin or tissue at the stoma site • Repetitive leakage of food or stomach contents • Pain, bleeding, puss, or inflammation at stoma site • Device is no longer fitting properly • Device falls out • Distended stomach

There is a risk that the device may be pulled from the stoma if undersized or excessive force was applied to device during use. There could be leakage of gastric content around the tube or leakage from the device. The tube could become clogged or have reduced flow. The interlock (feed port) may crack, leak, or separate from the device. The device may be difficult to remove. The device may become discolored over days to months of use.

Is the MiniONE® Non-Balloon Button MRI Compatible?

The MiniONE® Non-Balloon Button is considered MR Safe once placed in the stoma site when none of the MR Unsafe components are attached. The MiniONE® Non-Balloon Button, Placement T-Handle, feeding sets, syringes, and gauze are MR Safe. The Snap Arm Assembly Tool and Removal T-Handle are MR Unsafe.

How do you report adverse events?

- For Australian Users: Notice that any serious incident that occurs in relation to the device should be reported to Applied Medical Technology, Inc. and The Therapeutic Goods Administration (TGA) at <http://www.tga.gov.au>.
- Please contact AMT, our Authorized Representative (EC Rep), and/or the competent authority of the member state in which you are established if a serious incident has occurred in relation to the device.

Product Codes:

MiniONE® Non-Balloon Button Kits (Legacy)				MiniONE® Non-Balloon Button Kits (ENFit®)			
M1-1-1410	M1-1-1812	M1-1-2017	M1-1-2417	M1-1-1410-I	M1-1-1812-I	M1-1-2017-I	M1-1-2417-I
M1-1-1412	M1-1-1815	M1-1-2024	M1-1-2424	M1-1-1412-I	M1-1-1815-I	M1-1-2024-I	M1-1-2424-I
M1-1-1415	M1-1-1817	M1-1-2034	M1-1-2434	M1-1-1415-I	M1-1-1817-I	M1-1-2034-I	M1-1-2434-I
M1-1-1417	M1-1-1820	M1-1-2044	M1-1-2444	M1-1-1417-I	M1-1-1820-I	M1-1-2044-I	M1-1-2444-I
M1-1-1420	M1-1-1824			M1-1-1420-I	M1-1-1824-I		
M1-1-1425	M1-1-1830			M1-1-1425-I	M1-1-1830-I		
M1-1-1430	M1-1-1834			M1-1-1430-I	M1-1-1834-I		
	M1-1-1844				M1-1-1844-I		

Applied Medical Technology, Inc.

8006 Katherine Boulevard

Brecksville, OH 44141 USA

Toll Free: +1-800-869-7382 | Telephone: +1-440-717-4000

Website: www.appliedmedical.net | E-mail: cs@appliedmedical.net

© 2025 Applied Medical Technology, Inc. Rev C4746-C 07/2025

Patents: <http://www.AppliedMedical.net/Patents/> | Trademarks: <http://www.AppliedMedical.net/Trademarks/>

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