Patient Information Leaflet

AMT Traditional G-JET®

Traditional Length GJ Feeding Device

ENGLISH

What is in this leaflet?

This leaflet will answer some questions about the AMT Traditional G-JET®. AMT has provided this information as an educational resource tool. This is not intended as a substitute for professional medical care. Your FIRST source of information should be your healthcare provider.

What is the AMT Traditional G-JET®?

The AMT Traditional G-JET® provides for simultaneous gastric decompression / drainage and delivery of enteral nutrition into the distal duodenum or proximal jejunum. It enters the stomach through a gastric stoma. The tube is held in place (within the stoma tract) by means of an inflatable balloon and a sliding external bolster. The external tri-port contains two large ports; one labeled "JEJUNAL" and one labeled "GASTRIC". The JEJUNAL port is used for feeding into the small intestine. The GASTRIC port is used to drain the stomach with the use of low intermittent suction or gravity drainage. A third port labeled "BAL" is used to inflate and deflate the balloon. The device is available in various French sizes and lengths to accommodate the needs of patients.

The AMT Traditional G-JET® is made from medical-grade silicone (96%), medical-grade thermoplastic (2%), medical-grade silicone pad print ink (1%), and medical grade braided suture (1%). A stainless steel spring is contained within the jejunal tubing for the 16Fr and 18Fr devices. There are no manufacturing residuals that could pose a risk to the patient.

What is the AMT Traditional G-JET® Used for?

The AMT Traditional G-JET® is intended for gastric decompression / drainage and delivery of enteral nutrition into the distal duodenum or proximal jejunum. The device is indicated for use in adult, adolescent, child, and infant patients over 10kg who cannot absorb adequate nutrition through the stomach, who have intestinal motility problems, gastric outlet obstruction, severe gastroesophageal reflux, are at risk of aspiration, or in those who have had previous esophagectomy or gastrectomy. The use of this tube is also clinically indicated when simultaneous gastric decompression and jejunal feeding are needed. This includes patients in whom malnutrition already exists, or may result, secondary to concurrent conditions.

How do you use the AMT Traditional G-JET® after it has been placed?

Feeding: Open the cap to the jejunal access port (labeled "Jejunal") and connect the feed set to the tube. Deliver nutrition according to your healthcare professional's instructions. Make sure not to over-tighten the connection to the tube and ensure any clamps on the feed set are opened before starting to feed. If you see formula in the gastric drainage, stop feeding and notify your healthcare professional.

Flushing: Use room temperature water for tube flushing. The amount of water will depend on the patient's needs, clinical condition, and type of tube, but the average volume ranges from 10 to 50 ml for adults, and 3 to 10 ml for infants. Flush the feeding tube with water every 4-6 hours during continuous feeding, anytime the feeding is interrupted, before and after every intermittent feeding, or at least every 8 hours if the tube is not being used. Flush the feeding tube before and after medication channeling and between medications. Use a 30 to 60 ml syringe. Do not use smaller size syringes as this can increase pressure on the tube and potentially rupture smaller tubes. Do not use excessive force to flush the tube. Excessive force can perforate the tube and can cause injury to the gastrointestinal tract.

Medication Channeling: Use liquid medication when possible and consult the pharmacist to determine if it is safe to crush solid medication and mix with water. If safe, pulverize the solid medication into a fine powder form and dissolve the powder in water before channeling through the feeding tube. Never crush enteric coated medication or mix medication with formula. Using a syringe, flush the tube with the prescribed amount of water.


Figure 1: AMT Traditional G-JET®

Decompression: Some specialists recommend decompressing the stomach before or after feeding. Follow your specialist's instructions. The gastric port may be left open for gravity drainage or connected to low intermittent suction. Flush the gastric port EVERY SIX HOURS by injecting the prescribed amount of water. Do not use continuous or high intermittent suction. High pressure could collapse the tube or injure the stomach tissue and cause bleeding.

How do you take care of the AMT Traditional G-JET® after it has been placed?

Device Care: Check the device daily for any abnormalities such as damage, clogging or abnormal discoloration. Clogging and/or reduced flow are indications of diminished performance. The stoma area should be kept clean and dry at all times. It is important to clean the stoma site every day. A cotton swab or terry cloth can be used to clean the skin around the device with mild soap and water. Clean the stoma as directed by your healthcare professional. The feeding tube should be cleaned daily with warm water and mild soap, taking care not to pull or manipulate the tube excessively. The jejunal, gastric, and balloon ports should be cleaned daily using a cotton tip applicator or soft cloth to remove all residual formula and medication. DO NOT ROTATE THE TUBE. Jejunal tubes must not be rotated because they will kink and possibly lose position. Assess the patient for any signs of pain, pressure/discomfort, warmth, rashes, purulent or gastrointestinal drainage, pressure necrosis, skin breakdown, or hypergranulation tissue. If any of these symptoms are observed, contact your healthcare professional for advice.

Balloon Care: It is recommended that balloon volume is checked at least every two weeks, or as recommended by your healthcare professional. Remove the water with a syringe and compare the amount removed to the recommended amount. Refill the balloon and, if needed, add additional water to meet the recommended amount. Wait 10-20 minutes and repeat. The balloon is leaking if it has lost fluid, and the tube should be replaced. If the balloon is damaged, secure the tube in place using tape, then call your healthcare professional for instructions. Refill the balloon using sterile or distilled water, not air or saline. Saline can crystallize and clog the balloon valve or lumen, and air may seep out and cause the balloon to collapse. Be sure to use the recommended amount of water as over-inflation can obstruct the lumen or decrease balloon life and under-inflation will not secure the tube properly.

Clogging: First check to make sure that the tube is not kinked or clamped anywhere. If the clog is visible above the skin surface, attempt to massage the tube to break up the clog. Connect a 30 to 60 ml syringe filled with warm water into the appropriate adapter or lumen of the tube and gently push and pull the syringe plunger to free the clog. It may take several cycles of pushing / pulling the plunger to clear the clog. If clog cannot be removed, contact your healthcare professional, as the tube may need to be replaced. Do not use cranberry juice, cola drinks, meat tenderizer or chymotrypsin, as they can actually cause clogs or create adverse reactions in some patients.

How long does the AMT Traditional G-JET® last?

Balloon feeding devices are meant to be periodically replaced for optimal performance, functionality, and cleanliness. Device performance and functionality can degrade over time depending on usage and environmental conditions. Typical device longevity will vary for each patient depending on a number of factors, with typical device longevity ranging from 1-9 months. Some factors that can lead to reduced longevity include: gastric pH, diet of the patient, medications, balloon fill volume, trauma to the device, contact with sharp or abrasive objects, and overall tube care.

For optimal performance, it is recommended that the AMT Traditional G-JET® device be changed at least every 3 months or as often as indicated by your healthcare professional. Proactive replacement of the device will help ensure optimal functionality and will help prevent unexpected failure.

What are the possible side effects of having the AMT Traditional G-JET®?

The following complications may be associated with any transgastric-jejunal feeding device: Skin Breakdown • Infection • Hypergranulation Tissue • Stomach or Duodenal Ulcers • Intraperitoneal Leakage • Pressure Necrosis • Jejunal Limb Migration • Bowel perforation • Tubing Support (Spring) Dislodgement • Intussusception

Please consult your healthcare professional if you experience any of the following: Fever, vomiting, or diarrhea • Skin around the stoma site is red, discolored or raw • Drainage around the stoma site is white, yellow, or green; drainage has an unpleasant odor • Crusting is noted at the stoma site • Large amount of tissue build-up (such as granulation tissue) • Swollen skin or tissue at the stoma site • Repetitive leakage of food or stomach contents • Pain, bleeding, puss, or inflammation at stoma site • Device is no longer fitting properly • Device falls out • Distended stomach

There is a risk that if the internal balloon fails, the tube could fall out. There could be leakage of gastric content around the tube or leakage from the device. The tube could become clogged or have reduced flow. The device may become discolored over days to months of use.

Is the AMT Traditional G-JET® MRI Compatible?

Non-clinical testing has demonstrated that the AMT Traditional G-JET® feeding tube is MR Conditional. It can be scanned safely under the following conditions:

- Static magnetic field of 1.5-Tesla (1.5T) or 3.0-Tesla (3.0T).
- Spatial gradient field of up to:
 - 9,570 G/cm (95.70 T/m) for 1.5T systems
 - 5,720 G/cm (57.20 T/m) for 3.0T systems
- Maximum whole body averaged specific absorption rate (SAR) of:
 - 4.0 W/kg for 15 minutes of scanning in Normal Operating Mode at 1.5T.
 - 4.0 W/kg for 15 minutes of scanning in Normal Operating Mode at 3.0T.

1.5 RF heating: In non-clinical testing with body coil excitation, the AMT Traditional G-JET® feeding tube produced a temperature rise of less than 1.0C at a maximum whole body averaged specific absorption rate (SAR) of 4.0 W/kg, as assessed by calorimetry for 15 minutes of scanning in a 1.5T Siemens Espree (MRC30732) MR scanner with SYNGO MR B17 Software.

3.0 RF heating: In non-clinical testing with body coil excitation, the AMT Traditional G-JET® feeding tube produced a temperature rise of less than 1.0C at a maximum whole body averaged specific absorption rate (SAR) of 4.0 W/kg, as assessed by calorimetry for 15 minutes of scanning in a 3.0T Siemens Trio (MRC20587) MR scanner with SYNGO MR A30 4VA30A software.

3.0T MR Artifact: In testing using a 3.0T system with gradient-echo sequencing, the shape of the image artifact follows the approximate contour of the device and extends radially up to the 2.0 cm from the implant.

How do you report adverse events?

- For Australian Users: Notice that any serious incident that occurs in relation to the device should be reported to Applied Medical Technology, Inc. and The Therapeutic Goods Administration (TGA) at <http://www.tga.gov.au>.
- Please contact AMT, our Authorized Representative (EC Rep), and/or the competent authority of the member state in which you are established if a serious incident has occurred in relation to the device.

Product Codes:

AMT Traditional G-JET® Kits (Legacy)				AMT Traditional G-JET® Kits (ENFit®)			
14F	16F	18F	22F	14F	16F	18F	22F
TGJ-14-15	TGJ-16-15	TGJ-18-22	TGJ-22-30	TGJ-14-15-I	TGJ-16-15-I	TGJ-18-22-I	TGJ-22-30-I
TGJ-14-22	TGJ-16-22	TGJ-18-30	TGJ-22-45	TGJ-14-22-I	TGJ-16-22-I	TGJ-18-30-I	TGJ-22-45-I
TGJ-14-30	TGJ-16-30	TGJ-18-45		TGJ-14-30-I	TGJ-16-30-I	TGJ-18-45-I	
TGJ-14-45	TGJ-16-45			TGJ-14-45-I	TGJ-16-45-I		

Applied Medical Technology, Inc.

8006 Katherine Boulevard

Brecksville, OH 44141 USA

Toll Free: +1-800-869-7382 | Telephone: +1-440-717-4000

Website: www.appliedmedical.net | E-mail: cs@appliedmedical.net

© 2025 Applied Medical Technology, Inc. Rev C4748-C 07/2025

Patents: <http://www.AppliedMedical.net/Patents/> | Trademarks: <http://www.AppliedMedical.net/Trademarks/>

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