Patient Information Leaflet

AMT micro G-JET®

Micro Transgastric-Jejunal Feeding Device

ENGLISH

What is in this leaflet?

This leaflet will answer some questions about the AMT micro G-JET®. AMT has provided this information as an educational resource tool. This is not intended as a substitute for professional medical care. Your FIRST source of information should be your healthcare provider.

What is the AMT micro G-JET®?

The AMT micro G-JET® is a type of feeding tube that provides simultaneous gastric decompression / drainage and delivery of enteral nutrition into the small intestine (distal duodenum or proximal jejunum). It enters the stomach through an opening known as a stoma. The device includes a balloon that is filled with distilled or sterile water to help secure it inside the stomach and prevent it from falling out, along with an external bolster to help maintain tube position. The external bolster contains two ports; one labeled "JEJUNAL" (glow green) and one labeled "GASTRIC" (white). The JEJUNAL port is used for feeding into the small intestine while the GASTRIC port is used to drain the stomach with the use of low intermittent suction or gravity drainage. The micro G-JET® transitions from a 14F gastric segment to an 8F jejunal segment. The micro G-JET® is available in 14F with stoma lengths ranging from 0.8 cm to 2.5 cm with jejunal lengths ranging from 10 cm to 30 cm.

The AMT micro G-JET® is made from medical-grade silicone (64%), medical-grade thermoplastic (31%), stainless steel spring (4%), and medical-grade silicone pad print ink (1%). There are no manufacturing residuals that could pose a risk to the patient.

What is the AMT micro G-JET® Used for?

The AMT micro G-JET® is intended for gastric decompression / drainage and delivery of enteral nutrition into the distal duodenum or proximal jejunum. The device is indicated for use in adult, adolescent, child, and infant patients over 6 kg who cannot absorb adequate nutrition through the stomach, who have intestinal motility problems, gastric outlet obstruction, severe gastroesophageal reflux, are at risk of aspiration, or in those who have had previous esophagectomy or gastrectomy. The use of this tube is also clinically indicated when simultaneous gastric decompression and jejunal feeding are needed. This includes patients in whom malnutrition already exists, or may result, secondary to concurrent conditions.

How do you use the AMT micro G-JET® after it has been placed?

Feeding: Ensure the clamp is closed and attach the glow green adapter extension set into the glow green port labeled "JEJUNAL" by lining up the dark line on the extension set connector with the dark line on the jejunal feeding port. Lock the set into the jejunal feeding port by pushing in fully and rotating the connector CLOCKWISE until a slight resistance is felt (approximately 3/4 turn). Attach the opposite end of the extension set to the connector being used. Once connected, open the clamp to allow flow.

Figure 2: Attaching the Glow Green Extension Set

Medication Channeling: Use liquid medication when possible and consult the pharmacist to determine if it is safe to crush solid medication and mix with water. If safe, pulverize the solid medication into a fine powder form and dissolve the powder in water before channeling through the feeding tube. Never crush enteric coated medication or mix medication with formula. Using a syringe, flush the tube with the prescribed amount of water.


Figure 1: AMT micro G-JET®

Decompression: Some specialists recommend decompressing the stomach before or after feeding. Follow your specialist's instructions. Connect the AMT Gastric Extension Set to the GASTRIC port (white). The end of the Extension Set may be left open for gravity drainage or connected to low intermittent suction. Flush the gastric port EVERY FOUR HOURS by injecting the prescribed amount of water. Do not use continuous or high intermittent suction. High pressure could collapse the tube or injure the stomach tissue and cause bleeding.

Figure 3: Dual Extension/Feed Set Hook-Up

Flushing: Use room temperature water for tube flushing. The amount of water will depend on the patient's needs, clinical condition, and type of tube, but the average volume ranges from 3 to 10 ml for infants. Flush the feeding tube with water every 4 hours during continuous feeding, anytime the feeding is interrupted, before and after every intermittent feeding, or at least every 8 hours if the tube is not being used. Flush the feeding tube before and after medication channeling and between medications. A 5 ml syringe is recommended for volumes less than 5 ml. Caution should be taken if using a syringe smaller than 30 ml. If using a smaller syringe and resistance is felt, stop flushing immediately and switch to a larger syringe (30 to 60 ml). Continuing with a smaller syringe can increase pressure on the tube and potentially damage smaller tubes. Do not use excessive force to flush the tube. Excessive force can perforate the tube and can cause injury to the gastrointestinal tract.

How do you take care of the AMT micro G-JET® after it has been placed?

Device Care: Check the device daily for any abnormalities such as damage, clogging or abnormal discoloration. Clogging and/or reduced flow are indications of diminished performance. The stoma area should be kept clean and dry at all times. It is important to clean the stoma site every day. A cotton swab or terry cloth can be used to clean the skin around the device with mild soap and water. Clean the stoma as directed by your healthcare professional. The feeding tube should be cleaned daily with warm water and mild soap, taking care not to pull or manipulate the tube excessively. The jejunal, gastric, and balloon ports should be cleaned daily using a cotton tip applicator or soft cloth to remove all residual formula and medication.

Device Care: DO NOT ROTATE THE TUBE. Jejunal tubes must not be rotated because they will kink and possibly lose position. Assess the patient for any signs of pain, pressure/discomfort, warmth, rashes, purulent or gastrointestinal drainage, pressure necrosis, skin breakdown, or hypergranulation tissue. If any of these symptoms are observed, contact your healthcare professional for advice.

Balloon Care: It is recommended that balloon volume is checked at least every two weeks, or as recommended by your healthcare professional. Remove the water with a syringe and compare the amount removed to the recommended amount. Refill the balloon and, if needed, add additional water to meet the recommended amount. Wait 10-20 minutes and repeat. The balloon is leaking if it has lost fluid, and the tube should be replaced. If the balloon is damaged, secure the tube in place using tape, then call your healthcare professional for instructions.

Refill the balloon using sterile or distilled water, not air or saline. Saline can crystallize and clog the balloon valve or lumen, and air may seep out and cause the balloon to collapse. Be sure to use the recommended amount of water as over-inflation can obstruct the lumen or decrease balloon life and under-inflation will not secure the tube properly.

Clogging: First check to make sure that the tube is not kinked or clamped anywhere. If the clog is visible above the skin surface, attempt to massage the tube to break up the clog. Connect a 30 to 60 ml syringe filled with warm water into the appropriate adapter or lumen of the tube and gently push and pull the syringe plunger to free the clog. It may take several cycles of pushing / pulling the plunger to clear the clog. If clog cannot be removed, contact your healthcare professional, as the tube may need to be replaced. Do not use cranberry juice, cola drinks, meat tenderizer or chymotrypsin, as they can actually cause clogs or create adverse reactions in some patients.

How long does the AMT micro G-JET® last?

Low-profile balloon feeding devices are meant to be periodically replaced for optimal performance, functionality, and cleanliness. Device performance and functionality can degrade over time depending on usage and environmental conditions. Typical device longevity will vary for each patient depending on a number of factors, with typical device longevity ranging from 1-9 months. Some factors that can lead to reduced longevity include: gastric pH, diet of the patient, medications, balloon fill volume, trauma to the device, contact with sharp or abrasive objects, incorrect stoma length measurement, and overall tube care.

For optimal performance, it is recommended that the AMT micro G-JET® device be changed at least every 3 months or as often as indicated by your healthcare professional. Proactive replacement of the device will help ensure optimal functionality and will help prevent unexpected failure.

What are the possible side effects of having the AMT micro G-JET®?

The following complications may be associated with any transgastric-jejunal feeding device; any of which may result in potential patient injury, illness, and/or death. Skin Breakdown • Infection • Hypergranulation Tissue • Stomach or Duodenal Ulcers • Intraperitoneal Leakage • Pressure Necrosis • Intussusception • Gastric or Intestinal Perforation • Peritonitis • Device Migration

Please consult your healthcare professional if you experience any of the following: Fever, vomiting, or diarrhea • Skin around the stoma site is red, discolored or raw • Drainage around the stoma site is white, yellow, or green; drainage has an unpleasant odor • Crusting is noted at the stoma site • Large amount of tissue build-up (such as granulation tissue) • Swollen skin or tissue at the stoma site • Repetitive leakage of food or stomach contents • Pain, bleeding, puss, or inflammation at stoma site • Device is no longer fitting properly • Device falls out • Distended stomach

There is a risk that if the internal balloon fails, the tube could fall out. There could be leakage of gastric content around the tube or leakage from the device. The tube could become clogged or have reduced flow. The gastric or jejunal port may crack, leak, or separate from the device. The device may become discolored over days to months of use.

Is the AMT micro G-JET® MRI Compatible?

- Non-clinical testing has demonstrated that the AMT micro G-JET® is MR Conditional. It can be scanned safely under the following conditions:
 - Static magnetic field of 1.5-Tesla or 3-Tesla, only
 - Maximum spatial gradient magnetic field of 1,000-gauss/cm (10-T/m)
 - Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg in the Normal Operating Mode
- Under the scan conditions defined, the AMT micro G-JET® is expected to produce a maximum temperature rise of 2°C after 15 minutes of continuous scanning (i.e., per pulse sequence).
- In non-clinical testing, the image artifact caused by the AMT micro G-JET® extends approximately 20-mm from this device when imaged using a gradient echo pulse sequence and a 3-Tesla MR system.

How do you report adverse events?

- For Australian Users: Notice that any serious incident that occurs in relation to the device should be reported to Applied Medical Technology, Inc. and The Therapeutic Goods Administration (TGA) at <http://www.tga.gov.au>.
- Please contact AMT, our Authorized Representative (EC Rep), and/or the competent authority of the member state in which you are established if a serious incident has occurred in relation to the device.

Product Codes:

AMT micro G-JET® Kits (Legacy)			AMT micro G-JET® Kits (ENFit®)		
MGJ-1408-10	MGJ-1412-22	MGJ-1417-22	MGJ-1408-10-I	MGJ-1412-22-I	MGJ-1417-22-I
MGJ-1408-15	MGJ-1412-30	MGJ-1417-30	MGJ-1408-15-I	MGJ-1412-30-I	MGJ-1417-30-I
MGJ-1408-22	MGJ-1415-10	MGJ-1420-22	MGJ-1408-22-I	MGJ-1415-10-I	MGJ-1420-22-I
MGJ-1410-10	MGJ-1415-15	MGJ-1420-30	MGJ-1410-10-I	MGJ-1415-15-I	MGJ-1420-30-I
MGJ-1410-15	MGJ-1415-22	MGJ-1423-22	MGJ-1410-15-I	MGJ-1415-22-I	MGJ-1423-22-I
MGJ-1410-22	MGJ-1415-30	MGJ-1423-30	MGJ-1410-22-I	MGJ-1415-30-I	MGJ-1423-30-I
MGJ-1412-10	MGJ-1417-10	MGJ-1425-22	MGJ-1412-10-I	MGJ-1417-10-I	MGJ-1425-22-I
MGJ-1412-15	MGJ-1417-15	MGJ-1425-30	MGJ-1412-15-I	MGJ-1417-15-I	MGJ-1425-30-I

Applied Medical Technology, Inc.

8006 Katherine Boulevard, Brecksville, OH 44141 USA

Toll Free: +1-800-869-7382 | Telephone: +1-440-717-4000 | Website: www.appliedmedical.net | E-mail: cs@appliedmedical.net

© 2025 Applied Medical Technology, Inc. Rev C4749-C 07/2025

Patents: <http://www.AppliedMedical.net/Patents/> | Trademarks: <http://www.AppliedMedical.net/Trademarks/>

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