Patient Information Leaflet

AMT Balloon G-Tube

Replacement Gastrostomy Feeding Device

ENGLISH

What is in this leaflet?

This leaflet will answer some questions about the AMT Balloon G-Tube. AMT has provided this information as an educational resource tool. This is not intended as a substitute for professional medical care. Your FIRST source of information should be your healthcare provider.

What is the AMT Balloon G-Tube?

The AMT Balloon G-Tube is a gastrostomy feeding tube that is inserted directly into the stomach through an opening known as a stoma. The device includes a balloon that is filled with distilled or sterile water to help secure it inside the stomach and prevent it from falling out, along with a sliding external bolster to help maintain tube position. The sliding external bolster is adjustable to accommodate various stoma tract lengths. The AMT Balloon G-Tube is available in six sizes, 12F, 14F, 16F, 18F, 20F, and 24F, with either a Legacy or ENFit® connector.

The AMT Balloon G-Tube is made from medical-grade silicone (89%), medical-grade thermoplastic (10%), and medical-grade silicone pad print ink (1%). There are no manufacturing residuals that could pose a risk to the patient.

What is the AMT Balloon G-Tube Used for?

The AMT Balloon G-Tube is to be used as a percutaneous replacement gastrostomy tube for a patient with a well-established gastrostomy tract. This device assists in providing nutrition directly into the stomach through an established stoma in a human patient who is unable to consume nutrition by conventional means. The AMT Balloon G-Tube can also deliver medication and allow for decompression of the stomach.

How do you use the AMT Balloon G-Tube after it has been placed?

Feeding: Open the cap to the feed port and connect the feed set to the tube. Deliver nutrition according to your healthcare professional's instructions. Make sure not to over-tighten the connection to the tube and ensure any clamps on the feed set are opened before starting to feed.

Flushing: Use room temperature water for tube flushing. The amount of water will depend on the patient's needs, clinical condition, and type of tube, but the average volume ranges from 10 to 50 ml for adults, and 3 to 10 ml for infants. Flush the feeding tube with water every 4-6 hours during continuous feeding, anytime the feeding is interrupted, before and after every intermittent feeding, or at least every 8 hours if the tube is not being used. Flush the feeding tube before and after medication channeling and between medications. Use a 30 to 60 ml syringe. Do not use smaller size syringes as this can increase pressure on the tube and potentially rupture smaller tubes. Do not use excessive force to flush the tube. Excessive force can perforate the tube and can cause injury to the gastrointestinal tract.

Medication Channeling: Use liquid medication when possible and consult the pharmacist to determine if it is safe to crush solid medication and mix with water. If safe, pulverize the solid medication into a fine powder form and dissolve the powder in water before channeling through the feeding tube. Never crush enteric coated medication or mix medication with formula. Using a syringe, flush the tube with the prescribed amount of water.


Figure 1: AMT Balloon G-Tube

Decompression: Decompression can be performed if there is a need to check for residuals, or if the patient tends to retain air during feeding causing bloating and discomfort. To decompress, open the feed port and drain stomach contents into a container. After decompression is complete, flush with water.

How do you take care of the AMT Balloon G-Tube after it has been placed?

Device Care: Check the device daily for any signs of damage or clogging. Clogging and/or reduced flow are indications of diminished performance. The stoma area should be kept clean and dry at all times. It is important to clean the stoma site every day. A cotton swab or terry cloth can be used to clean the skin around the device with mild soap and water. Clean the stoma as directed by your healthcare professional. The device should be rotated daily for site hygiene. Always check the stoma site for redness, pain/soreness, swelling, or any drainage. If any of these symptoms are observed, contact your healthcare professional for advice.

Balloon Care: It is recommended that the balloon volume is checked at least every two weeks, or as instructed by your healthcare professional. Remove the water with a syringe and compare the amount removed to the recommended amount. Refill the balloon and, if needed, add additional water to meet the recommended amount. Wait 10-20 minutes and repeat. The balloon is leaking if it has lost fluid, and the tube should be replaced. If the balloon is damaged, secure the tube in place using tape, then call your healthcare professional for instructions.

Refill the balloon using sterile or distilled water, not air or saline. Saline can crystallize and clog the balloon valve or lumen, and air may seep out and cause the balloon to collapse. Be sure to use the recommended amount of water as over-inflation can obstruct the lumen or decrease balloon life and under-inflation will not secure the tube properly.

Clogging: First check to make sure that the tube is not kinked or clamped anywhere. If there is a visible clog in the tubing, attempt to massage the device to break up the clog. Connect a 30 to 60 ml syringe filled with warm water into the appropriate adapter or lumen of the tube and gently push and pull the syringe plunger to free the clog. It may take several cycles of pushing / pulling the plunger to clear the clog. If clog cannot be removed, contact your healthcare professional, as the tube may need to be replaced. Do not use excessive force or pressure to attempt to clear the clog. This can cause the tubing to rupture.

How long does the AMT Balloon G-Tube last?

Balloon feeding devices are meant to be periodically replaced for optimal performance, functionality, and cleanliness. Device performance and functionality can degrade over time depending on usage and environmental conditions. Typical device longevity will vary for each patient depending on a number of factors, with typical device longevity ranging from 1-9 months. Some factors that can lead to reduced longevity include: gastric pH, diet of the patient, medications, balloon fill volume, trauma to the device, contact with sharp or abrasive objects, and overall tube care.

For optimal performance, it is recommended that the AMT Balloon G-Tube device be changed at least every 3 months or as often as indicated by your healthcare professional. Proactive replacement of the device will help ensure optimal functionality and will help prevent unexpected failure.

What are the possible side effects of having the AMT Balloon G-Tube?

Potential complications when using the AMT Balloon G-Tube include but are not limited to: Aspiration • Abscess, wound infection and skin breakdown • Hypergranulation tissue • Buried bumper syndrome • Pressure necrosis • Gastrointestinal bleeding and/or ulcerations • Ileus or gastroparesis • Intraperitoneal leakage • Bowel and gastric volvulus • Peritonitis • Gastrocolic fistula • Sepsis • Obstruction

Please consult your healthcare professional if you experience any of the following: Nausea, vomiting, abdominal bloating, or diarrhea • Pain, bleeding, and/or inflammation at the G-tube site • Crusting at the stoma site • Skin around the stoma site that is red, discolored, or raw • Stoma site drainage and/or puss that is white, yellow, or green and may smell bad • Repetitive leakage of food or stomach contents • Distended stomach • Fever • Tube clog • Migration of the tubing resulting in the device extending longer or shorter from the stoma • Extubation resulting in the device becoming removed from the stoma and unable to be replaced easily • Distinct indentation at the G-tube site or a distinct gap between the device and skin

There is a risk that if the internal balloon fails, the tube could fall out. There could be leakage of gastric content around the tube or leakage from the device. The tube could become clogged or have reduced flow. The device may become discolored over days to months of use.

Is the AMT Balloon G-Tube MRI Compatible?

The AMT Balloon G-Tube is considered MR Safe once placed.

How do you report adverse events?

- For Australian Users: Notice that any serious incident that occurs in relation to the device should be reported to Applied Medical Technology, Inc. and The Therapeutic Goods Administration (TGA) at <http://www.tga.gov.au>.
- Please contact AMT, our Authorized Representative (EC Rep), and/or the competent authority of the member state in which you are established if a serious incident has occurred in relation to the device.

Product Codes:

AMT Balloon G-Tube (Legacy)	AMT Balloon G-Tube (ENFit®)
7-1210	7-1210-I
7-1410	7-1410-I
7-1610	7-1610-I
7-1810	7-1810-I
7-2010	7-2010-I
7-2410	7-2410-I

Applied Medical Technology, Inc.
8006 Katherine Boulevard
Brecksville, OH 44141 USA

Toll Free: +1-800-869-7382 | Telephone: +1-440-717-4000
Website: www.appliedmedical.net | E-mail: cs@appliedmedical.net

© 2025 Applied Medical Technology, Inc. Rev C4750-C 07/2025

Patents: <http://www.AppliedMedical.net/Patents/> | Trademarks: <http://www.AppliedMedical.net/Trademarks/>

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