# RHAPSO®

## Suture Passer

- ES** RHAPSO®. Pasador de sutura.
- DE** RHAPSO®. Fadentransporter-Nadel.
- IT** RHAPSO®. Passasutura.
- FR** RHAPSO®. Passe-suture.
- NO** RHAPSO®. Suture Passer.
- SV** RHAPSO®. Suturförare.
- PT-BR** RHAPSO®. Passador de Sutura.
- NL** RHAPSO®. Hechtdraadgrijper.
- DA** RHAPSO®. Suturpasser.
- FI** RHAPSO®. Ommelohjain.
- AR** RHAPSO®. ممسك غرز
- BG** RHAPSO®. Уред за прокараване на шевове.
- CS** RHAPSO®. Podávač sutury.
- ET** RHAPSO®. Suture Passer.
- EL** RHAPSO®. Συσκευή Διέλευσης Ραμμάτων.
- HU** RHAPSO®. Varróanyag-átvezető.
- JA** RHAPSO®. 縫合糸パッサー.
- LV** RHAPSO®. Šuvju padevējs.
- LT** RHAPSO®. Siūlo įvediklis.
- PT** RHAPSO®. Passador de sutura.
- PL** RHAPSO®. Przeciagacz do nici.
- RO** RHAPSO®. Instrument de trecere a suturii.
- SK** RHAPSO®. Podávač sutúry.
- SL** RHAPSO®. Priprava za podajanje materiala za šivanje.
- TR** RHAPSO®. Sütür Geçirici.


STERILE EO

DEHP

MD Rx ONLY

<p><b>ENGLISH</b></p> <p><b>Figure 1: RHAPSO®</b></p> <ol style="list-style-type: none"> <li>1. Straight Needle</li> <li>2. Graspers</li> <li>3. Curved Needle</li> <li>4. Instrument Body</li> <li>5. Luer Hub Cap</li> <li>6. Luer Hub</li> </ol>	<p><b>ESPAÑOL</b></p> <p><b>Figura 1: RHAPSO®</b></p> <ol style="list-style-type: none"> <li>1. Aguja recta</li> <li>2. Pinzas</li> <li>3. Aguja curvada</li> <li>4. Cuerpo del instrumento</li> <li>5. Tapón del conector Luer</li> <li>6. Conector Luer</li> </ol>	<p><b>DEUTSCH</b></p> <p><b>Abbildung 1: RHAPSO®</b></p> <ol style="list-style-type: none"> <li>1. Gerade Nadel</li> <li>2. Greifer</li> <li>3. Gebogene Nadel</li> <li>4. Instrumentenkörper</li> <li>5. Luer-Hülsenverschluss</li> <li>6. Luer-Hülse</li> </ol>
<p><b>ITALIANO</b></p> <p><b>Figura 1: RHAPSO®</b></p> <ol style="list-style-type: none"> <li>1. Ago diritto</li> <li>2. Pinze</li> <li>3. Ago curvo</li> <li>4. Corpo dello strumento</li> <li>5. Tappo dell'attacco Luer</li> <li>6. Attacco Luer</li> </ol>	<p><b>FRANÇAIS</b></p> <p><b>Figure 1: RHAPSO®</b></p> <ol style="list-style-type: none"> <li>1. Aiguille droite</li> <li>2. Pinces</li> <li>3. Aiguille incurvée</li> <li>4. Corps de l'instrument</li> <li>5. Bouchon obturateur pour raccord conique Luer</li> <li>6. Raccord Luer</li> </ol>	<p><b>NORSK</b></p> <p><b>Figur 1: RHAPSO®</b></p> <ol style="list-style-type: none"> <li>1. Rett nål</li> <li>2. Gripere</li> <li>3. Kurvet nål</li> <li>4. Instrumentets hovedenhet</li> <li>5. Lokk til luer hub</li> <li>6. Luer hub</li> </ol>
<p><b>SVENSKA</b></p> <p><b>Bild 1: RHAPSO®</b></p> <ol style="list-style-type: none"> <li>1. Rak nål</li> <li>2. Gripare</li> <li>3. Böjd nål</li> <li>4. Instrumentkropp</li> <li>5. Skydd till luerhätta</li> <li>6. Luerhätta</li> </ol>	<p><b>PORTUGUÊS</b></p> <p><b>Figura 1: RHAPSO®</b></p> <ol style="list-style-type: none"> <li>1. Agulha reta</li> <li>2. Pinças</li> <li>3. Agulha curvada</li> <li>4. Corpo do instrumento</li> <li>5. Tampa do Luer Hub</li> <li>6. Luer Hub</li> </ol>	<p><b>NEDERLANDS</b></p> <p><b>Afbeelding 1: RHAPSO®</b></p> <ol style="list-style-type: none"> <li>1. Rechte naald</li> <li>2. Grippers</li> <li>3. Gebogen naald</li> <li>4. Behuizing van instrument</li> <li>5. Dop van luer-hub</li> <li>6. Luer-hub</li> </ol>
<p><b>DANSK</b></p> <p><b>Figur 1: RHAPSO®</b></p> <ol style="list-style-type: none"> <li>1. Lige kanyle</li> <li>2. Gribetænger</li> <li>3. Krum kanyle</li> <li>4. Instrumenthovedel</li> <li>5. Hætte til luer-nav</li> <li>6. Luer-nav</li> </ol>	<p><b>SUOMI</b></p> <p><b>Kuva 1: RHAPSO®</b></p> <ol style="list-style-type: none"> <li>1. Suora neula</li> <li>2. Tartuntapihdit</li> <li>3. Kaareva neula</li> <li>4. Instrumentin runko</li> <li>5. Luer Hub -korkki</li> <li>6. Luer Hub</li> </ol>	<p>عربي</p> <p><b>الشكل 1: RHAPSO®</b></p> <ol style="list-style-type: none"> <li>1. إبرة مستقيمة</li> <li>2. ممسكات</li> <li>3. إبرة منحنية</li> <li>4. جسم الأداة</li> <li>5. غطاء وصلات لور</li> <li>6. وصلة لور</li> </ol>

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## Instructions for Use

**Caution:** Federal (USA) law restricts this device to sale by or on the order of a physician. This device is provided STERILE for single use. Please inspect all contents for damage. If damaged or sterile barrier is breached, do not use.

### DEVICE DESCRIPTION

The RHAPSO® is a sterile, single-case use, hand-held suture grasping device that is designed to pass suture through soft tissue. The device consists of an instrument body with an attached hypodermic needle and a slidable actuated grasper that is housed within the lumen of the needle.

### INTENDED USE

The RHAPSO® is intended to pass suture through soft tissues of the body. This product is intended for use by surgeons trained in endoscopic and/or laparoscopic techniques and/or interventional radiologists.

**WARNING: THIS DEVICE IS INTENDED FOR SINGLE USE. DO NOT REUSE, REPROCESS OR RE-STERILIZE THE MEDICAL DEVICE. DOING SO MAY COMPROMISE BIOCOMPATIBILITY CHARACTERISTICS, DEVICE PERFORMANCE AND/OR MATERIAL INTEGRITY; ANY OF WHICH MAY RESULT IN POTENTIAL PATIENT INJURY, ILLNESS AND/OR DEATH.**

### INDICATIONS FOR USE

The RHAPSO® is indicated to pass suture through soft tissues of the body during endoscopic/laparoscopic surgery and interventional radiology procedures in child, adolescent and adult populations.

### CONTRAINDICATIONS FOR USE

There are no contraindications for use of the RHAPSO®.

### COMPLICATIONS

At this time, there are no complications associated with the use of the RHAPSO®.

**NOTE:** Please contact AMT, our Authorized Representative (EC Rep), and/or the competent authority of the member state in which you are established if a serious incident has occurred in relation to the device.

### CLINICAL BENEFITS, PERFORMANCE CHARACTERISTICS

Clinical benefits to be expected when using the RHAPSO® include but are not limited to:

- Small needle (17G) is designed to minimize trauma and enhance wound healing
- Luer-lock compatible hub with patent lumen enables delivery and removal of fluid or radiographic media
- Unique grasper arm reduces the complexity of intracorporeal suture placement, enhancing procedural efficiency
- Small needle (17G) with straight and curved options increases flexibility in suture placement and needle approach

Performance characteristics of the RHAPSO® include but are not limited to:

- Device provides the ability to pass and retrieve sutures through soft tissues
- Device can securely grasp a broad range of suture sizes (5-0 to 0)
- Straight and curved echogenic needle options for procedural versatility
- Ergonomic design provides controlled, reliable performance for seamless surgical precision

### INSTRUCTIONS

**WARNING: THIS DEVICE IS TO BE USED BY TRAINED HEALTHCARE PROFESSIONALS AND SHOULD NOT BE USED CLINICALLY WITHOUT A FULL UNDERSTANDING OF THE DIRECTIONS FOR USE.**

**WARNING: INJURY TO TISSUE CAN OCCUR FROM CONTACT WITH THE POINTED END OF THE NEEDLE.**


**NOTE:** Safe use of the device requires proper visualization of the working end of the needle. Non-clinical testing has shown that the needle and grasper can be sufficiently visualized under endoscopic, laparoscopic, fluoroscopic, or ultrasound guidance.

**NOTE:** Safe use of the device requires an adequate intracorporeal space to deploy the grasper to the open position. The grasper must be allowed to fully deploy in order to easily grasp suture.

**NOTE:** The grasper components of the device may become damaged if overloaded with force. To avoid damage to the grasper components, do not exceed the tensile strength of the suture.

**WARNING: RESORBABLE MONOFILAMENT SUTURES SUCH AS CATGUT, MAY BE DAMAGED WHEN CAPTURED BY THE GRASPER. IT IS RECOMMENDED THAT THESE SUTURES BE CAPTURED WITH THE DEVICE AND INSPECTED FOR DAMAGE PRIOR TO BEING USED IN THE PATIENT.**

The grasper of the RHAPSO® has two operating positions – open and closed. The grasper is opened by depressing the luer hub. The grasper is closed automatically by an internal return spring.

## INSTRUCTIONS

### 1. Introduce the RHAPSO® and deliver the suture.

The grasper should be in the closed position with the suture captured before introducing the RHAPSO®. With the grasper in the open position, envelop the suture such that the grasper mechanism ensnares the suture. To capture the suture and close the grasper, release the luer hub. The device should return to the closed position automatically and the suture should be captured.

**NOTE:** Cut and remove needle from suture, if present. Ensure only the suture strand is to be passed.



With the grasper in the closed position, use the RHAPSO® with captured suture to percutaneously access the intracorporeal site where suture is to be delivered.

**NOTE:** A small stab incision through the skin may aid in needle insertion.

Verify internal position of the RHAPSO® before proceeding.

**NOTE:** Contrast medium can be injected through the RHAPSO® to aid in confirming position under fluoroscopy.



**NOTE:** Luer access may also be used for aspiration or fluid delivery (channeling) through the device.

### 2. Release the Suture and Withdraw the RHAPSO®.

Open the grasper to release the suture by depressing the luer hub. Maneuver the RHAPSO® away from the suture so that it is not immediately re-captured. The internal spring will return the grasper to the closed position.

The grasper is now in the closed position and may be withdrawn from the patient, leaving the suture through the tissue.


### 3. Reintroduce the RHAPSO® and Capture the Suture.

Reintroduce the RHAPSO® while still in the closed position.

**NOTE:** If the suture knot is to be buried subcutaneously, the RHAPSO® should be re-introduced through the same skin-level incision as the suture and directed at an angle to achieve the desired suture bite width.

Once in the lumen, open the grasper and maneuver the grasper to envelop the suture such that the grasper mechanism ensnares the suture.

Release the luer hub so that the grasper returns to the closed position, securing the suture.


### 4. Pass the Suture.

With the suture captured, carefully withdraw the RHAPSO® from the patient; the suture is pulled in tow.

**NOTE:** With the grasper closed, the suture is mechanically trapped inside the needle. However, if excessive tension is applied asymmetrically to one of the suture limbs, the suture may slide out of the grasper.

Once the instrument is clear of the patient's tissue, open the grasper to release the suture.


**WARNING: THE MECHANISM OF THE GRASPER MAY RESULT IN PERMANENT DEFORMATION OF MONOFILAMENT SUTURES (KINKING) THAT MAY COMPROMISE THE TENSILE STRENGTH OF THE SUTURE. FOR SAFE USE OF THE DEVICE, IT IS RECOMMENDED THAT ANY PORTION OF THE SUTURE THAT IS DEFORMED BE REMOVED FROM INSIDE THE PATIENT.**

### 5. Secure the Suture.

Tension and knot the suture. Cut and discard the excess suture from the ends.

**NOTE:** The device can be disposed of by following facility protocol and in accordance with local, state, and federal laws and regulations.

## MRI SAFETY INFORMATION

The RHAPSO® is MR-Unsafe. The RHAPSO® is not to be used in or around the MR environment. Suture placed using the RHAPSO® is MR-safe. A patient with suture may be scanned safely.

## THANK YOU!

Thank you for choosing AMT. For additional help and information regarding the use of our device, feel free to contact AMT with the contact information on the back of the instructions for use. We are happy to hear your thoughts and help with your concerns and questions.

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Applied Medical Technology, Inc.