




## AMT G-JET® Low Profile Gastric-Jejunal Enteral Tube

### Product Sizes and HCPC Codes

#### AMT G-JET® Balloon Buttons Size Chart 18F Coming Soon!

14 French			16 French		
Order #	Stoma Length cm	Jejunal Length cm	Order #	Stoma Length cm	Jejunal Length cm
GJ-1410-15	1.0	15	GJ-1610-15	1.0	15
GJ-1410-22	1.0	22	GJ-1610-22	1.0	22
GJ-1412-15	1.2	15	GJ-1612-15	1.2	15
GJ-1412-22	1.2	22	GJ-1612-22	1.2	22
GJ-1412-30	1.2	30	GJ-1612-30	1.2	30
GJ-1415-15	1.5	15	GJ-1615-15	1.5	15
GJ-1415-22	1.5	22	GJ-1615-22	1.5	22
GJ-1415-30	1.5	30	GJ-1615-30	1.5	30
GJ-1415-45	1.5	45	GJ-1615-45	1.5	45
GJ-1417-15	1.7	15	GJ-1617-15	1.7	15
GJ-1417-22	1.7	22	GJ-1617-22	1.7	22
GJ-1417-30	1.7	30	GJ-1617-30	1.7	30
GJ-1417-45	1.7	45	GJ-1617-45	1.7	45
GJ-1420-22	2.0	22	GJ-1620-22	2.0	22
GJ-1420-30	2.0	30	GJ-1620-30	2.0	30
GJ-1420-45	2.0	45	GJ-1620-45	2.0	45
GJ-1423-30	2.3	30	GJ-1623-30	2.3	30
GJ-1423-45	2.3	45	GJ-1623-45	2.3	45
GJ-1425-30	2.5	30	GJ-1625-30	2.5	30
GJ-1425-45	2.5	45	GJ-1625-45	2.5	45
GJ-1427-30	2.7	30	GJ-1627-30	2.7	30
GJ-1427-45	2.7	45	GJ-1627-45	2.7	45
GJ-1430-45	3.0	45	GJ-1630-45	3.0	45
GJ-1435-45	3.5	45	GJ-1635-45	3.5	45



**F = Tube Diameter**  
**cm = Gastric Length**  
**cm = Jejunal Length**

Example: **14F** x **2.3cm** x **30cm**

(F) French Size  
 (cm) Centimeter

Available in 14F  
 The smallest GJ in the world  
 Coming Soon - 18F G-JET®

\*Included with device:  
 Jejunal (feed set) 8-1257G 12" RA Glow Green connector with improved Y-Port  
 Gastric (extension set) 6-1222 12" RA white connector with clear Y-Port  
 All G-JET Devices are packaged 1/Box

NOTE: G-JET devices should never be cut to custom lengths

#### AMT G-JET® Balloon Buttons Fill Volume Table

French Size	Minimum	Recommended	Maximum
14F	3 ml	4 ml	5 ml
16F	3 ml	4 ml	5 ml

#### HCPC Codes for Enteral Feeding Device and Accessories

##### Caregiver/ Patient Ordering Information

AMT takes orders directly but does not process medical billing paperwork. Prepaid/credit card orders are accepted with a prescription.

French Size	Description	HCPC Codes	Description
M1-X-XXXX	MiniONE® Button Low Profile Feeding Device Family	B4088	Low Profile G-Tube
GJ-XXXX-XX	G-JET® Button Low Profile Gastric-Jejunal Enteral Tube	B4088	Low Profile GJ-Tube
8-XXXX	MiniONE® Button Feed Set Extensions	B9998	Generic Accessory Feed Set
5-XXXX	MINI Button Low Profile Balloon Feeding Device	B4088	Low-Profile G-Tube
6-XXXX	MINI Button Feed Set Extensions	B9998	Generic Accessory Feed Set
7-XXXX	AMT Balloon Replacement Gastrostomy Tube	B4087	Standard G-Tube

AMT has provided the information above as a guide only and makes no guarantee as to the accuracy and reliability of the information: it is not a comprehensive resource. We strongly suggest contacting your local Medicaid or Medicare office to determine correct billing procedures, reimbursement fees, and policies. AMT assumes no responsibility for errors, omissions, or changes in reimbursement policies.

\*The AMT G-JET® has mutually exclusive ports which require two different feed/extension sets for improved patient safety:

- The gastric extension set will only fit the gastric (white) port.
- The jejunal feed set will only fit the jejunal (Glow Green) port.

Please Call or Email Customer Service at:  
 Phone- (800 869 7382) Email- CS@AppliedMedical.net

Manufactured by: Applied Medical Technology, Inc.



Applied Medical Technology, Inc.  
 8006 Katherine Boulevard  
 Brecksville, Oh 44141 USA  
 Phone: 800 869 7382  
 Fax: 440 717 4200  
 Email: CS@AppliedMedical.net  
 www.AppliedMedical.net  
 ISO 13485 Registered Company



Rev 4 - 08/15