



## AMT MiniONE® Balloon Button

### Product Sizes and HCPC Codes

#### AMT MiniONE® Balloon Buttons

12 French		14 French		16 French		18 French		20 French		24 French	
Order #	Stoma Length (cm)	Order #	Stoma Length (cm)	Order #	Stoma Length (cm)	Order #	Stoma Length (cm)	Order #	Stoma Length (cm)	Order #	Stoma Length (cm)
M1-5-1208	0.8	M1-5-1408	0.8								
M1-5-1210	1.0	M1-5-1410	1.0	M1-5-1610	1.0	M1-5-1810	1.0	M1-5-2010	1.0		
M1-5-1212	1.2	M1-5-1412	1.2	M1-5-1612	1.2	M1-5-1812	1.2	M1-5-2012	1.2		
M1-5-1215	1.5	M1-5-1415	1.5	M1-5-1615	1.5	M1-5-1815	1.5	M1-5-2015	1.5	M1-5-2415	1.5
M1-5-1217	1.7	M1-5-1417	1.7	M1-5-1617	1.7	M1-5-1817	1.7	M1-5-2017	1.7	M1-5-2417	1.7
M1-5-1220	2.0	M1-5-1420	2.0	M1-5-1620	2.0	M1-5-1820	2.0	M1-5-2020	2.0	M1-5-2420	2.0
M1-5-1223	2.3	M1-5-1423	2.3	M1-5-1623	2.3	M1-5-1823	2.3	M1-5-2023	2.3	M1-5-2423	2.3
M1-5-1225	2.5	M1-5-1425	2.5	M1-5-1625	2.5	M1-5-1825	2.5	M1-5-2025	2.5	M1-5-2425	2.5
M1-5-1227	2.7	M1-5-1427	2.7	M1-5-1627	2.7	M1-5-1827	2.7	M1-5-2027	2.7	M1-5-2427	2.7
M1-5-1230	3.0	M1-5-1430	3.0	M1-5-1630	3.0	M1-5-1830	3.0	M1-5-2030	3.0	M1-5-2430	3.0
M1-5-1235	3.5	M1-5-1435	3.5	M1-5-1635	3.5	M1-5-1835	3.5	M1-5-2035	3.5	M1-5-2435	3.5
M1-5-1240	4.0	M1-5-1440	4.0	M1-5-1640	4.0	M1-5-1840	4.0	M1-5-2040	4.0	M1-5-2440	4.0
		M1-5-1444	4.4	M1-5-1644	4.4	M1-5-1844	4.4				
		M1-5-1450	5.0								
		M1-5-1455	5.5			M1-5-1855	5.5				
		M1-5-1460	6.0								
		M1-5-1465	6.5								

#### AMT MiniONE® Balloon Button Fill Volume Table

French Size	Minimum	Recommended	Maximum
12F	2 ml	2.5 ml	3 ml
14F	3 ml	4 ml	5 ml
16F	4 ml	6 ml	8 ml
18F	6 ml	8 ml	10 ml
20F	7 ml	10 ml	15 ml
24F	7 ml	10 ml	15 ml


#### HCPC Codes for Enteral Feeding Devices

##### Caregiver/ Patient Ordering Information

AMT takes orders directly but does not process medical billing paperwork. Prepaid /credit card orders are accepted with a prescription.

French Size	Description	HCPC Codes	Description
M1-X-XXXX	MiniONE® Button Low Profile Feeding Device Family	B4088	Low Profile G-Tube
GJ-XXXX-XX	G-JET® Button Low Profile Gastric-Jejunal Enteral Tube	B4088	Low Profile GJ-Tube
8-XXXX	MiniONE® Button Feed Set Extensions	B9998	Generic Accessory Feed Set
5-XXXX	MINI Button Low Profile Balloon Feeding Device	B4088	Low-Profile G-Tube
6-XXXX	MINI Button Feed Set Extensions	B9998	Generic Accessory Feed Set
7-XXXX	AMT Balloon Replacement Gastrostomy Tube	B4087	Standard G-Tube

AMT has provided the information above as a guide only and makes no guarantee as to the accuracy and reliability of the information: it is not a comprehensive resource. We strongly suggest contacting your local Medicaid or Medicare office to determine correct billing procedures, reimbursement fees, and policies. AMT assumes no responsibility for errors, omissions, or changes in reimbursement policies.



F = Tube Diameter  
 cm = Stoma Length

Example: 14F x 2.3 cm

(F) French Size  
 (cm) Centimeter



Applied Medical Technology, Inc.  
 8006 Katherine Boulevard  
 Brecksville, Oh 44141 USA  
 Phone: 800 869 7382  
 Fax: 440 717 4200  
 Email: CS@AppliedMedical.net  
 www.AppliedMedical.net  
 ISO 13485 Registered Company



Please Call or Email Customer Service at:  
 Phone- (800 869 7382) Email- CS@AppliedMedical.net