

AMT Bridle™

Nasal Tube Retaining System Family



Reduce Pullouts + Decrease Cost

- Bridling is more effective than tape¹.
- Bridling reduced additional KUB x-rays/radiographs by 99.3%².

Improve Nutrition + Recovery

- Bridled patients have seen >23 uninterrupted feeding days³.
- Consistent nutrition delivery can minimize length-of-stay⁴.

Patient Safety + Comfort

- New Probe design has a 20% smaller diameter⁵.
- Pre-Attached Clip for increased safety and improved handling.
- Monofilament Tubing is easy to clean and comfortable for patients.

Pro Range Clips

- Clips accommodate several nasal tube sizes.



5F-6F




8F-10F



12F-14F



16F-18F




Innovating. Educating. Changing Lives.™

AMT Offers the Largest Range of Pediatric + Adult Sizes⁵

Bridle Pro® and Pro Range Systems

An AMT Exclusive. The New Standard of Care.

Bridle Pro® Range clips accommodate several nasal tube sizes which means fewer products to keep in stock.



Bridle Pro® Range Clips:



AMT Bridle™ & Bridle Pro® Clips:




AMT Bridle™ Safe + Secure

The AMT Bridle™ features a superior clip: it is textured for an easier grip, it is pre-attached for improved handling, and it has rounded edges for less irritation.

AMT Bridle™ has a smaller probe than alternative brands for a more comfortable placement.

Large Bore Nasal Tube Securement:
14F-18F can secure Salem Sump® tubes.



AMT has the ONLY Bridle systems FDA cleared for pediatric patients.

Pediatric Bridle Pro® Range		
Description	Order #	Box/Qty
5F-6F	4-420506P	5
8F-10F	4-420810P	5

Bridle Pro® Range: Blue Bridle Tubing

Pediatric Bridle Pro®		
Description	Order #	Box/Qty
8F	4-4208P	5

Bridle Pro®: Blue Bridle Tubing

Bridle Pro® Range		
Description	Order #	Box/Qty
8F-10F	4-420810	5
12F-14F	4-421214	5
16F-18F	4-421618	5

Bridle Pro® Range: Blue Bridle Tubing

Bridle Pro®		
Description	Order #	Box/Qty
8F	4-4208	5
10F	4-4210	5
12F	4-4212	5
14F*	4-4214	5
16F*	4-4216	5
18F*	4-4218	5

Bridle Pro®: Blue Bridle Tubing

AMT Bridle™		
Description	Order #	Box/Qty
8F	4-4108	5
10F	4-4110	5
12F	4-4112	5
14F*	4-4114	5
16F*	4-4116	5
18F*	4-4118	5

AMT Bridle™: White Umbilical Tape

Kit Contents: Flexible Bridle Catheter, Retrieval Probe, Pre-Attached Clip, Safety Stylet Guide, Removal Tool, Packet of Water Soluble Lubricant.

References: [1] Bechtold, M.L., et al. Nutr Clin Pract. 2014;29:667-671. [2] Gupta, P.K., et al. Gastroenterology 2010;138(5):S-235. [3] Power, S., et al., e-SPEN 2010;5:e73-e76. [4] Badosa, E.L. Nutr Hosp. 2017;34(4):907-913. [5] Based on internal testing.

*14F-18F can be used with a Salem Sump®

For a complete list of products visit us at www.AppliedMedical.net
800 869 7382 | CS@AppliedMedical.net | ICS@AppliedMedical.net



AMT Bridle™ Family



Applied Medical Technology, Inc.
8006 Katherine Boulevard
Brecksville, OH 44141
www.AppliedMedical.net
An ISO 13485 Registered Company
MR Safe (with stylet removed)
Made In USA



© 2023 Applied Medical Technology, Inc. C-3802-P

Patents: <http://www.AppliedMedical.net/Patents/> | Trademarks: <http://www.AppliedMedical.net/Trademarks/>
Salem Sump® is a registered trademark of Cardinal Health Ireland Unlimited Company.