

AMT MaxVac™

Surgical Suction Instrument



Universal Compatibility

Designed to connect to most O.R. suction systems with no capital investment.

Multiple Tip Options

Easily exchangeable - choose the vented tip to fit your needs, any time.

Intuitive Assembly + Use

Device is ready to use in seconds and features an ambidextrous design.

Clears Fluids with Ease

Suctions liquids of different viscosities quickly and effectively.

Large Volume Body Fluid Extractor



Innovating. Educating. Changing Lives.™

AMT MaxVac™



Surgical Suction Instrument

High Viscosity Body Fluid Removal


The AMT MaxVac™ is a unique tool that quickly removes body fluids during surgery. With multiple tip options, the adaptable device can be used for a variety of surgical procedures.

Device Anatomy


- 1 Inlet Tip - 8 mm Diameter
- 2 Inlet Tip - 18 mm Diameter
- 3 Nozzle (28 mm Diameter)
- 4 Reservoir (500 ml Capacity)
- 5 Thumb Throttle + Plug
- 6 Alignment Indicator
- 7 No-Slip Design
- 8 O.R. Suction Adapter




Assembly + Use*




Connect




Assemble



Attach



Insert + Extract



Tilt + Drain

Possible Applications:

- Blood (Hematoma)
- Pus
- Amniotic Fluid
- Bowel Contents
- Seroma
- Vomitus
- Peritoneal Fluid
- Other Body Fluids

AMT MaxVac™

ORDER NUMBER	BOX/QTY
ST7002	1

Packed sterile, with 2 tip options.

For a complete list of products visit us at www.AppliedMedical.net
800 869 7382 | CS@AppliedMedical.net | ICS@AppliedMedical.net



AppliedMedical.net



Applied Medical Technology, Inc.
8006 Katherine Boulevard
Brecksville, OH 44141
www.AppliedMedical.net
An ISO 13485 Registered Company
Made In USA

