

RHAPSO™

Versatile Suture Passing Instrument

NEW

Small Gauge Suture Passer

17-gauge needle designed for minimally invasive procedures.

Streamlined + Secure


Wide opening easily captures a broad range of sutures.

Preferred Approach

Fine needle allows great accessibility to tight spaces.

Patent Lumen

Luer access allows delivery and removal of fluid or radiographic media.



Precision in Every Pass



Innovating. Educating. Changing Lives.™

800 869 7382 | CS@AppliedMedical.net | ICS@AppliedMedical.net


RHAPSO™

Versatile Suture Passing Instrument

Simplified Suturing


Engineered for precision and agility, the intuitive design of **RHAPSO™** enables confident suture passing and retrieving. Its spring-assisted retriever utilizes atraumatic grasper arms, ideal for holding sutures securely. The device's ergonomic features provide controlled, reliable performance for seamless surgical outcomes.

- 1 - Luer Accessible Hub Grasper Activator
- 2 - Device Housing
- 3 - Needle (17 Ga)
Straight or Curved Options, Echogenic
- 4 - Grasper Arms




RHAPSO™ | Order Information

STRAIGHT NEEDLE		CURVED NEEDLE	
SD110-S-07-101	BOX 1	SD110-C-07-101	BOX 1
SD110-S-07-105	BOX 5	SD110-C-07-105	BOX 5




Align Opening

*Needle images are transparent for clarity only.



Precise Grasp



Secure Suture

For a complete list of products, visit us at www.AppliedMedical.net
800 869 7382 | CS@AppliedMedical.net | ICS@AppliedMedical.net



www.AppliedMedical.net



Applied Medical Technology, Inc.
 8006 Katherine Boulevard
 Brecksville, OH 44141
www.AppliedMedical.net
 An ISO 13485 Registered Company
 Made In USA

