


GEMINI[®]

Suture Passer with Magnet Technology

NEW



Self-Aligning Magnet Technology

Easily navigate suture placement within the surgical field.

Preferred Suture Compatibility

Loop feature allows for exchange to suture of choice.

Patent Lumen

Luer access allows delivery and removal of fluid and radiographic media.

Smallest Gauge Suture Passer

17-gauge needle designed to minimize trauma and enhance wound healing.

Elegance, Ease, & Versatility



Innovating. Educating. Changing Lives.™

800 869 7382 | CS@AppliedMedical.net | ICS@AppliedMedical.net

GEMINI®

Suture Passer with Magnet Technology

Make Connections


The GEMINI® is designed to pass suture through soft tissue during minimally invasive procedures. The magnetized arm of the grasper works in conjunction with a magnetic suture to reduce the technical difficulty of intracorporeal suture retrieval.

GEMINI® | Order Information


ORDER NUMBER	QTY/BOX
CURVED NEEDLE KITS*	
SD221-C-07-101	1
SD221-C-07-105	5
STRAIGHT NEEDLE KITS*	
SD221-S-07-101	1
SD221-S-07-105	5
SUTURE LOOPS	
S4-8000	10

*Kit contents: GEMINI® device, introducer needle, two magnetic suture loops


- 1 - Luer Accessible Hub
- 2 - Activator Slide
- 3 - Suture (3-0, Braided Polyester) in Introducer Tube (2)
- 4 - Needle (17 Ga) and Type Straight or Curved
- 5 - Retractable End Effector Grasper Arm, Magnet, Trap




Prompt Placement




PLACE



CONNECT



TRAP



RETRIEVE

For a complete list of products, visit us at www.AppliedMedical.net
800 869 7382 | CS@AppliedMedical.net | ICS@AppliedMedical.net



AppliedMedical.net



Applied Medical Technology, Inc.
8006 Katherine Boulevard
Brecksville, OH 44141
www.AppliedMedical.net
An ISO 13485 Registered Company
Made In USA



© 2025 Applied Medical Technology, Inc. C-4996-D

Patents: <http://www.AppliedMedical.net/Patents/> | Trademarks: <http://www.AppliedMedical.net/Trademarks/>

*This document is not a substitute for the IFU/DFU (Instructions for Use/Directions for Use).