


# Comprehensive I.P. Kits

## Initial Placement Gastrostomy Sets

The intuitive design of AMT I.P. Kits ensures all of the essential tools for initial placement of a percutaneous gastrostomy feeding device are packaged in one easy-to-use kit.



### Quality, Right from the Start

#### 1 Introducer Needles

Two lengths (short, long) to accommodate a broad range of patients.

#### 2 Introducer / Stiffener

An AMT exclusive, size-specific, guidewire compatible.

#### 3 Slip Tip Syringe

For use with Balloon Stoma Measuring Device.

#### 4 BSMD

Guidewire compatible Balloon Stoma Measuring Device.

#### 5 Safety Scalpel

#### 6 Gauze

#### 7 Guidewire

With J-tip

#### 8 Multiple Size Dilators

### Designed to Minimize Trauma Associated with Stoma Tract Development

- Tapered Dilators Help Ensure Proper Dilation
- Requires Only One Endoscopic Pass
- Compatible with Laparoscopic, Endoscopic, and Fluoroscopic Placement
- For Use with 10F - 18F Devices



Innovating. Educating. Changing Lives.™

# Comprehensive I.P. Kits

Intuitive Design. AMT Quality.

AMT offers Comprehensive and Basic I.P. Kits that are designed to place 10F - 18F gastrostomy devices. Comprehensive kits include dilators in a range of French sizes, a guidewire, two introducer needles (short, long), a balloon stoma measuring device, an introducer/stiffener, a scalpel, a small syringe, and gauze. Basic kits include dilators in a range of French sizes, a guidewire, two introducer needles (short, long), and an introducer/stiffener.

## Comprehensive I.P. Kits<sup>‡</sup>

ORDER #	KIT CONTENTS
IP-10-S	(3) Dilators (7F, 12F, <b>14F-S</b> )
	(1) Guidewire, 70 cm (0.018")
	(2) Introducer Needles, 18 Gauge (short, long)
	(1) Guidewire Compatible BSMD
	(1) Guidewire Compatible Introducer/Stiffener
	(1) Safety Scalpel, #11
	(1) 5 cc Syringe
	(4) 4" x 4" Gauze
ORDER #	KIT CONTENTS
IP-1218	(5) Dilators (7F, 12F, 16F, 18F, 22F)
	(1) Guidewire, 90 cm (0.038")
	(2) Introducer Needles, 18 Gauge (short, long)
	(1) Guidewire Compatible BSMD
	(2) Guidewire Compatible Introducer/Stiffener
	(1) Safety Scalpel, #11
	(1) 5 cc Syringe
	(4) 4" x 4" Gauze
ORDER #	KIT CONTENTS
IP-1218-S	(5) Dilators (7F, 12F, 16F, <b>18F-S</b> , 22F)
	(1) Guidewire, 90 cm (0.038")
	(2) Introducer Needles, 18 Gauge (short, long)
	(1) Guidewire Compatible BSMD
	(2) Guidewire Compatible Introducer/Stiffener
	(1) Safety Scalpel, #11
	(1) 5 cc Syringe
	(4) 4" x 4" Gauze

I.P. Kits are packaged 1/box (sterile). | <sup>‡</sup>NOTE: Kits do not contain T-Fasteners or U-Stitches.  
"-S" I.P. Kits include one sheathed dilator, noted in bold above.







## Basic I.P. Kits<sup>\*</sup>

ORDER #	KIT CONTENTS
IP-DIL-10	(3) Dilators (7F, 12F, 14F)
	(1) Guidewire, 45 cm (0.018")
	(2) Introducer Needles, 18 Gauge (short, long)
	(1) Guidewire Compatible Introducer/Stiffener
ORDER #	KIT CONTENTS
IP-DIL	(4) Dilators (7F, 12F, 16F, 20F)
	(1) Guidewire, 45 cm (0.038")
	(2) Introducer Needles, 18 Gauge (short, long)
	(1) Guidewire Compatible Introducer/Stiffener

I.P. Kits are packaged 1/box (sterile). | \* NOTE: Kits do not contain T-Fasteners or U-Stitches.

## I.P. Kit Accessories

Available Separately

ORDER #	DESCRIPTION
1804	Balloon Stoma Measuring Device Guidewire Compatible (up to 0.038") 
1704	Balloon Stoma Measuring Device Non-Guidewire Compatible 
4-6010	Introducer/Stiffener Guidewire Compatible (up to 0.018") Device Compatibility: 10F and larger Low Profile Devices 
4-6000	Introducer/Stiffener - 4.6" Length Guidewire Compatible (up to 0.038") Device Compatibility: 12F and larger Low Profile Devices 
4-6200	Introducer/Stiffener - 9.6" Length Guidewire Compatible (up to 0.038") Device Compatibility: 12F and larger Low Profile Devices 
4-5000	Introducer/Stiffener Non-Guidewire Compatible Device Compatibility: 12F and larger MINI™ Button Devices 

I.P. Kit accessories are individually pouched and sold 10/box (sterile).

For a complete list of products visit us at [www.AppliedMedical.net](http://www.AppliedMedical.net)  
800 869 7382 | [CS@AppliedMedical.net](mailto:CS@AppliedMedical.net) | [ICS@AppliedMedical.net](mailto:ICS@AppliedMedical.net)



Comprehensive I.P. Kits



Applied Medical Technology, Inc.  
8006 Katherine Boulevard  
Brecksville, OH 44141  
[www.AppliedMedical.net](http://www.AppliedMedical.net)  
An ISO 13485 Registered Company  
Cleared For Initial Placement  
Made In USA

