

Proper Fit **Every Time**

To ensure a proper fit, it is crucial to measure the stoma length using the AMT Balloon Stoma Measuring Device.

Balloon Stoma Measuring Device	Order #	Qty/Box
Guidewire Compatible	1804	10
Non-Guidewire Compatible	1704	10

G-JET® Button

Low Profile Gastric-Jejunal Enteral Feeding Tube

G-JET® Button

G	J	14 F	16 F	18 F
LENGTH		LEGACY ORDER NUMBERS*		
1.0 cm	15 cm	GJ-1410-15	GJ-1610-15	
	22 cm	GJ-1410-22	GJ-1610-22	
1.2 cm	15 cm	GJ-1412-15	GJ-1612-15	GJ-1812-22
	22 cm	GJ-1412-22	GJ-1612-22	
	30 cm	GJ-1412-30	GJ-1612-30	GJ-1812-30
1.5 cm	15 cm	GJ-1415-15	GJ-1615-15	
	22 cm	GJ-1415-22	GJ-1615-22	GJ-1815-22
	30 cm	GJ-1415-30	GJ-1615-30	GJ-1815-30
	45 cm	GJ-1415-45	GJ-1615-45	GJ-1815-45
1.7 cm	15 cm	GJ-1417-15	GJ-1617-15	
	22 cm	GJ-1417-22	GJ-1617-22	GJ-1817-22
	30 cm	GJ-1417-30	GJ-1617-30	GJ-1817-30
2.0 cm	45 cm	GJ-1417-45	GJ-1617-45	GJ-1817-45
	22 cm	GJ-1420-22	GJ-1620-22	GJ-1820-22
	30 cm	GJ-1420-30	GJ-1620-30	GJ-1820-30
2.3 cm	45 cm	GJ-1420-45	GJ-1620-45	GJ-1820-45
	22 cm	GJ-1423-22		
	30 cm	GJ-1423-30	GJ-1623-30	GJ-1823-30
2.5 cm	45 cm	GJ-1423-45	GJ-1623-45	GJ-1823-45
	22 cm	GJ-1425-22		
	30 cm	GJ-1425-30	GJ-1625-30	GJ-1825-30
2.7 cm	45 cm	GJ-1425-45	GJ-1625-45	GJ-1825-45
	30 cm	GJ-1427-30	GJ-1627-30	GJ-1827-30
	45 cm	GJ-1427-45	GJ-1627-45	GJ-1827-45
3.0 cm	45 cm	GJ-1430-45	GJ-1630-45	GJ-1830-45
3.5 cm	45 cm	GJ-1435-45	GJ-1635-45	GJ-1835-45
4.0 cm	45 cm	GJ-1440-45	GJ-1640-45	GJ-1840-45
4.4 cm	45 cm	GJ-1444-45	GJ-1644-45	GJ-1844-45
5.0 cm	45 cm	GJ-1450-45	GJ-1650-45	GJ-1850-45
5.5 cm	45 cm	GJ-1455-45	GJ-1655-45	GJ-1855-45
6.0 cm	45 cm	GJ-1460-45	GJ-1660-45	GJ-1860-45
6.5 cm	45 cm	GJ-1465-45	GJ-1665-45	GJ-1865-45
7.0 cm	45 cm	GJ-1470-45	GJ-1670-45	GJ-1870-45


The G-JET® Button is sold 1/box (sterile).
NOTE: G-JET® devices should never be cut to custom lengths.
Guidewire available (180 cm): Order Number 4-GW-5
*Add "-I" for ENFit® G-JET® Button kits.

Balloon Fill Volumes

The recommended water volume is clearly printed on the fill valve of each G-JET® Button. Consult your healthcare provider before adjusting the fill volume.


FRENCH SIZE	MIN	REC.	MAX
14F	3 ml	4 ml	5 ml
16F	3 ml	4 ml	5 ml
18F	7 ml	8 ml	9 ml

Extension + Feed Sets Included in Kit



ADAPTERS	RIGHT ANGLE CONNECTORS			STRAIGHT CONNECTORS		
	2"	12"	24"	2"	12"	24"
COMPATIBLE GASTRIC EXTENSION SETS						
Bolus		6-1221	6-2421	6-0211	6-1211	6-2411
Clear Y-Port		6-1222	6-2422	6-0212	6-1212	
ENFit® Y-Port		6-1222-ISOSAF				
COMPATIBLE JEJUNAL FEED SETS						
Purple Y-Port		8-1257G	8-2455			
ENFit® Y-Port		8-1255	8-2455	8-0212		
		8-1255-ISOSAF	8-2455-ISOSAF	8-0212-ISOSAF		
Bolus		8-1221	8-2421	8-0211	8-1211	8-2411
ENFit® Single Port		8-1221-ISOSAF	8-2421-ISOSAF	8-0211-ISOSAF	8-1211-ISOSAF	8-2411-ISOSAF
Clear Y-Port XL		8-0222	8-2422	8-0212	8-1212	

Legacy Included ENFit® Included Extension + Feed sets come packaged as 10/box.



Reduce Clogging

Anti-Kink Technology

AMT Exclusive Balloon

Improve Comfort



Applied Medical Technology, Inc.
8006 Katherine Boulevard
Brecksville, OH 44141
www.AppliedMedical.net
An ISO 13485 Registered Company
Cleared For Initial Placement
MR Conditional
Made In USA




© 2025 Applied Medical Technology, Inc. C-4497-N
Patents: <http://www.AppliedMedical.net/Patents/> | Trademarks: <http://www.AppliedMedical.net/Trademarks/>
*The Anti-Reflux Valve is opened when the feed/extension set is locked into place. | ENFit® is a registered trademark of GEDSA.
1.Claim based on internal testing of 14F-22F devices. | 2.Claim based on internal testing of 16F devices.

AMT logo

Innovating. Educating. Changing Lives.™

Same French Size, Larger Lumens

DESIGNED TO REDUCE CLOGGING




An Advantage in Enteral Nutrition

Anti-Kink. Exclusive AMT Balloon. Glow Green™ Technology.

Reduce Leaks, Superior Comfort

- 1 Anti-Reflux Valve**
Prevents backflow from the G + J ports²
 - 2 Tapered Transition**
For a more secure fit
 - 3 Exclusive AMT Balloon**
Apple-shaped, expands horizontally
Creates a superior gastric seal²
- Soft + Flexible**
Medical grade silicone is gentle to the skin and easy to clean around
- Four Point Stabilization**
Reduces motion that can cause granulation tissue



Safety + Innovation

- 4 Guidewire Adapter**
For ease of placement
- 5 Mutually Exclusive G+J**
Lowers risk of misconnecting feed and extension sets
- 6 Glow Green™ J-Port**
Easier nighttime feeds
Feed the Green!
- 7 Dual Safety Straps**
Allow independent access to the G and/or J port
- 8 Anti-Kink Technology**
16F and 18F sizes offer spring reinforcement to maintain flow

SUCCESS STORIES

My G-JET® Lasts Longer

The length of time between routine exchanges has significantly increased since I switched to the G-JET®. Interventional radiology is now replacing my G-JET® every 4-5 months... because it can last!

Amanda's G-JET®
AppliedMedical.net/Amandas-G-JET

Mutually Exclusive Ports Leading the Way in Patient Safety

Accidentally feeding the stomach instead of the jejunum can be dangerous! Mutually exclusive **GASTRIC** and **JEJUNAL** ports are designed to reduce the risk of misconnections. Different feed/extension sets are required for each port, so it is imperative to maintain adequate stock of both sets. Always remember to **Feed the Green!**

