

# Explant Express®

Breast Implant Suction Retrieval Device

**NEW**



## Rapid Removal

Avoids spillage of ruptured implants – keeping silicone gel away from gloves, instruments, and surgical field.

## Intuitive Design

Effortless performance requires no advanced training. Ergonomic no-slip design is suitable for all hand sizes.

## Cost Effective

Universal compatibility with standard operating room suction. Rapid retrieval reduces operating time.

**Avoid the mess with Explant Express®!**



Innovating. Educating. Changing Lives.™

800 869 7382 | [CS@AppliedMedical.net](mailto:CS@AppliedMedical.net) | [ICS@AppliedMedical.net](mailto:ICS@AppliedMedical.net)

# Explant Express<sup>®</sup>

## Breast Implant Suction Retrieval Device

### Engineered for Simplicity

Redefining the standard of care, the Explant Express<sup>®</sup> streamlines the breast implant removal process. This cost-effective solution features a universal design, suitable for a variety of breast implant sizes. While traditional removal methods often take up to an hour, the Explant Express<sup>®</sup> can remove a ruptured implant in a matter of seconds!<sup>1</sup> The device's separable halves allow for post-procedural inspection of extracted material.



**MESS!**



**EXPRESS!**

Avoid the mess with Explant Express<sup>®</sup>!

Product Number: **ST7001**  
Packaged 1/Box (Sterile).

For a complete list of products, visit us at [www.AppliedMedical.net](http://www.AppliedMedical.net)  
800 869 7382 | [CS@AppliedMedical.net](mailto:CS@AppliedMedical.net) | [ICS@AppliedMedical.net](mailto:ICS@AppliedMedical.net)



Explant Express<sup>®</sup>



Applied Medical Technology, Inc.  
8006 Katherine Boulevard  
Brecksville, OH 44141  
[www.AppliedMedical.net](http://www.AppliedMedical.net)  
An ISO 13485 Registered Company  
Made In USA



© 2026 Applied Medical Technology, Inc. **C-5108-A**

Patents: <http://www.AppliedMedical.net/Patents/> | Trademarks: <http://www.AppliedMedical.net/Trademarks/>

\*This document is not a substitute for the IFU/DFU (Instructions for Use/Directions for Use). | 1. Based on internal testing.