

















AMT PRODUCTS

MR Status Quick Reference Guide

Always consult the product-specific Instructions for Use (IFU) and institutional MRI safety protocols prior to scanning. Final responsibility for patient safety rests with the MR operator and supervising clinician.

ENTERAL PRODUCTS:

NutraGlide®		MR Safe (once stylet is removed)
AMT Bridle Pro®		MR Safe (once probe, catheter, and stylet is removed)
micro G-JET®		MR Conditional - consult IFU for scan conditions
G-JET® Button		MR Conditional - consult IFU for scan conditions
Traditional G-JET®		MR Conditional - consult IFU for scan conditions
MiniONE® Balloon Button		MR Conditional - consult IFU for scan conditions
MiniONE® Capsule Button		MR Safe
MiniONE® Non-Balloon Button		MR Safe
Traditional Balloon G-Tube		MR Safe
Traditional Capsule Monarch® G-Tube		MR Safe
Traditional Capsule Dome G-Tube		MR Safe
Decompression Tube		MR Safe
MiniACE® Button		MR Conditional - consult IFU for scan conditions
MiniSPC® Button		MR Conditional - consult IFU for scan conditions
MINI™ Classic Button		MR Conditional - consult IFU for scan conditions
CINCH®		MR Safe

SURGICAL PRODUCTS:

IntelAir®		MR Safe
-----------	---	---------

For a complete list of products, visit us at www.AppliedMedical.net
800 869 7382 | CS@AppliedMedical.net | ICS@AppliedMedical.net



AMT



Applied Medical Technology, Inc.
8006 Katherine Boulevard
Brecksville, OH 44141
An ISO 13485 Registered Company
Made In USA

