

AMT Bridle Pro[®]

Nasal Tube Retaining System



Reduce Pullouts + Decrease Cost

- Bridling is more effective than tape.¹
- Bridling reduced additional KUB x-rays/radiographs by 99.3%.²

Improve Nutrition + Recovery

- Bridled patients have seen >23 uninterrupted feeding days.³
- Consistent nutrition delivery can minimize length-of-stay.⁴

Patient Safety + Comfort

- New probe design has a 20% smaller diameter.⁵
- Pre-attached clip for increased safety and improved handling.
- Monofilament tubing is easy to clean and comfortable for patients.

Pro Range Clips

- Clips accommodate several nasal tube sizes.



5F - 6F




8F - 10F



12F - 14F



16F - 18F




Innovating. Educating. Changing Lives.™

800 869 7382 | CS@AppliedMedical.net | ICS@AppliedMedical.net

AMT Offers the Largest Range of Pediatric + Adult Sizes⁵

AMT Bridle Pro® and Pro Range Systems An AMT Exclusive. The New Standard of Care.

AMT Bridle Pro® Range clips accommodate several nasal tube sizes, which means fewer products to keep in stock.



AMT Bridle Pro® Safe + Secure

The AMT Bridle Pro® features a superior clip; it is textured for an easier grip, pre-attached for improved handling, and has rounded edges for less irritation.

AMT Bridle Pro® has a smaller probe than alternative brands for a more comfortable placement.

AMT Bridle Pro® Range Clips:



AMT Bridle Pro® Clips:



Large Bore Nasal Tube Securement:

14F-18F can secure Salem Sump® tubes.



AMT has the ONLY Bridle system FDA cleared for pediatric patients.

AMT Bridle Pro® | Order Information

Size	Pediatric AMT Bridle Pro® Range	Pediatric AMT Bridle Pro®	AMT Bridle Pro® Range	AMT Bridle Pro®	Qty/Box
5F	4-420506P				5
6F					5
8F	4-420810P	4-4208P	4-420810	4-4208	5
10F				4-4210	5
12F			4-421214	4-4212	5
14F*				4-4214	5
16F*			4-421618	4-4216	5
18F*				4-4218	5

Kit Contents: Flexible bridle catheter, retrieval probe, pre-attached clip, safety stylet guide, removal tool, packet of water-soluble lubricant.

References: [1] Bechtold, M.L., et al. Nutr Clin Pract. 2014;29:667-671. [2] Gupta, P.K., et al. Gastroenterology 2010;138(5):S-234-S-235. [3] Power, S., et al., e-SPEN 2010;5:e73-e76. [4] Badosa, E.L. Nutr Hosp. 2017;34(4):907-913. [5] Based on internal testing.

*14F-18F can be used with a Salem Sump®.

For a complete list of products, visit us at www.AppliedMedical.net
800 869 7382 | CS@AppliedMedical.net | ICS@AppliedMedical.net



AMT Bridle Pro®



Applied Medical Technology, Inc.
8006 Katherine Boulevard
Brecksville, OH 44141
www.AppliedMedical.net
An ISO 13485 Registered Company
Made In USA



© 2026 Applied Medical Technology, Inc. C-3802-T

Patents: <http://www.AppliedMedical.net/Patents/> | Trademarks: <http://www.AppliedMedical.net/Trademarks/>
Salem Sump® is a registered trademark of Cardinal Health Ireland Unlimited Company.