

AMT
MiniACE® Button
Cecostomy + Appendicostomy

MiniACE® Button

Low Profile Antegrade Continence Enema Button

MiniACE® Enema Button | Order Information

Length (cm)	ORDER NUMBERS					
	10F		12F		14F	
	LEGACY	ENFIT®	LEGACY	ENFIT®	LEGACY	ENFIT®
0.8	ACE-1008	ACE-1008-I	ACE-1208	ACE-1208-I	ACE-1408	ACE-1408-I
1.0	ACE-1010	ACE-1010-I	ACE-1210	ACE-1210-I	ACE-1410	ACE-1410-I
1.2	ACE-1012	ACE-1012-I	ACE-1212	ACE-1212-I	ACE-1412	ACE-1412-I
1.5	ACE-1015	ACE-1015-I	ACE-1215	ACE-1215-I	ACE-1415	ACE-1415-I
1.7	ACE-1017	ACE-1017-I	ACE-1217	ACE-1217-I	ACE-1417	ACE-1417-I
2.0	ACE-1020	ACE-1020-I	ACE-1220	ACE-1220-I	ACE-1420	ACE-1420-I
2.3	ACE-1023	ACE-1023-I	ACE-1223	ACE-1223-I	ACE-1423	ACE-1423-I
2.5	ACE-1025	ACE-1025-I	ACE-1225	ACE-1225-I	ACE-1425	ACE-1425-I
2.7	ACE-1027	ACE-1027-I	ACE-1227	ACE-1227-I	ACE-1427	ACE-1427-I
3.0	ACE-1030	ACE-1030-I	ACE-1230	ACE-1230-I	ACE-1430	ACE-1430-I
3.5	ACE-1035	ACE-1035-I	ACE-1235	ACE-1235-I	ACE-1435	ACE-1435-I
4.0	ACE-1040	ACE-1040-I	ACE-1240	ACE-1240-I	ACE-1440	ACE-1440-I
4.4	ACE-1044	ACE-1044-I	ACE-1244	ACE-1244-I	ACE-1444	ACE-1444-I
5.0	ACE-1050	ACE-1050-I	ACE-1250	ACE-1250-I	ACE-1450	ACE-1450-I
5.5	ACE-1055	ACE-1055-I	ACE-1255	ACE-1255-I	ACE-1455	ACE-1455-I
6.0	ACE-1060	ACE-1060-I	ACE-1260	ACE-1260-I	ACE-1460	ACE-1460-I
6.5	ACE-1065	ACE-1065-I	ACE-1265	ACE-1265-I	ACE-1465	ACE-1465-I
7.0	ACE-1070	ACE-1070-I	ACE-1270	ACE-1270-I	ACE-1470	ACE-1470-I
8.0	ACE-1080	ACE-1080-I	ACE-1280	ACE-1280-I	ACE-1480	ACE-1480-I
9.0	ACE-1090	ACE-1090-I	ACE-1290	ACE-1290-I	ACE-1490	ACE-1490-I
10.0	ACE-10100	ACE-10100-I	ACE-12100	ACE-12100-I	ACE-14100	ACE-14100-I

MiniACE® Button comes packaged as Qty/Box 1 (sterile). | Kit contents include: device (1), irrigation set (1), 5cc syringe (1), introducer/stiffener (1), gauze, and water soluble lubricant (1). MR Conditional.

Irrigation Set Included in Kit



Complete Portfolio – Irrigation Sets

ADAPTER OPTIONS	RIGHT ANGLE CONNECTOR OPTIONS	
	12"	24"
Bolus	2-1221	2-2421
	2-1221-ISOSAF	2-2421-ISOSAF

Note: Irrigation sets are packaged as Qty/Box 10 (non-sterile)



MiniACE® Button



© 2026 Applied Medical Technology, Inc. C-4642-G

Patents: <http://www.AppliedMedical.net/Patents/> | Trademarks: <http://www.AppliedMedical.net/Trademarks/>

Chait Trapdoor™ is a trademark of Cook® Medical | ENFit® is a registered trademark of GEDSA. | 1.Claim based on internal testing of 10F-14F devices.

Balloon Fill Volumes

FRENCH SIZE	MIN	REC.	MAX
10F	1 ml	1.5 ml	2 ml
12F	2 ml	2.5 ml	3 ml
14F	3 ml	4 ml	5 ml

Recommended water fill volume is clearly printed on the balloon fill valve of each MiniACE® Button. Consult your healthcare provider before adjusting fill volumes.


Balloon Measuring Device | Order Number 1904

Guidewire Compatible (up to 0.038"), 10F Tubing




AMT CINCH® Tube Securement

Order Numbers: CINCH410S, CINCH417M, and CINCH930L



AMT Clamp™ | Order Number 4-3000



AppliedMedical.net/Stories/Audreys-MiniACE-Enema-Button/



Cecostomy + Appendicostomy

Cleared for Initial Placement

Superior Leakage Control

Improved Stoma Retention

Applied Medical Technology, Inc.
8006 Katherine Boulevard
Brecksville, OH 44141
www.AppliedMedical.net
An ISO 13485 Registered Company
Cleared For Initial Placement
Made In USA



Innovating. Educating. Changing Lives.™

Secure Attachment



Size Variety




Comfortable




MiniACE® Button

- Cleared for Initial Placement
- Minimizes stoma dilation associated with conversion from Chait Trapdoor™ to MiniACE®
- Cecostomy + Appendicostomy



Cecostomy




Appendicostomy

Patient Comfort First

- Medical grade silicone, gentle to the skin
- Flexible + soft, easy to clean around the button
- Exclusive AMT Balloon, designed to reduce pullouts and leakage
- Balloon bolster shields the internal tip of the MiniACE®, protecting the cecum


Reduce Leaks, Superior Retention

- AMT's MiniACE® Button is 450% more secure within the stoma when compared to the Chait Trapdoor™.¹
- Exclusive AMT Balloon, along with tapered tubing at the external bolster, reduce the risk of leaks.
- Anti-leakage valve helps prevent backflow.




Safe and Secure Irrigation Set Attachment


- The connection between the MiniACE® Button and its irrigation set is 10x stronger than the competition.¹




STABILIZE



ALIGN & INSERT



ROTATE



LOCK



IRRIGATE