

MiniONE® Button

Capsule Non-Balloon and Non-Balloon Devices



Exclusive Capsule Technology

Material facilitates easy and comfortable insertion.

Glow Green™ Technology

Easier nighttime feeds.

Improve Comfort

Combines the advantage of a low profile device with the longer life of a non-balloon device.

Non-Balloon Advantage

Ideal for patients who don't tolerate balloon buttons.

More Durable

Internal bolster is smaller than most devices, durable and less likely to block the pylorus.¹



Pull for Capsule Deployment



Innovating. Educating. Changing Lives.™


800 869 7382 | CS@AppliedMedical.net | ICS@AppliedMedical.net

MiniONE® Capsule Non-Balloon & Non-Balloon Devices

MiniONE® | Order Information

	LENGTH (cm)	14F		18F		20F		24F	
		LEGACY ORDER #	ENFIT® ORDER #	LEGACY ORDER #	ENFIT® ORDER #	LEGACY ORDER #	ENFIT® ORDER #	LEGACY ORDER #	ENFIT® ORDER #
Non-Balloon Button	1.0	M1-1-1410	M1-1-1410-I						
	1.2	M1-1-1412	M1-1-1412-I	M1-1-1812	M1-1-1812-I				
	1.5	M1-1-1415	M1-1-1415-I	M1-1-1815	M1-1-1815-I				
	1.7	M1-1-1417	M1-1-1417-I	M1-1-1817	M1-1-1817-I	M1-1-2017	M1-1-2017-I	M1-1-2417	M1-1-2417-I
	2.0	M1-1-1420	M1-1-1420-I	M1-1-1820	M1-1-1820-I				
	2.4			M1-1-1824	M1-1-1824-I	M1-1-2024	M1-1-2024-I	M1-1-2424	M1-1-2424-I
	2.5	M1-1-1425	M1-1-1425-I						
	3.0	M1-1-1430	M1-1-1430-I	M1-1-1830	M1-1-1830-I				
	3.4			M1-1-1834	M1-1-1834-I	M1-1-2034	M1-1-2034-I	M1-1-2434	M1-1-2434-I
Capsule Non-Balloon Button	4.4			M1-1-1844	M1-1-1844-I	M1-1-2044	M1-1-2044-I	M1-1-2444	M1-1-2444-I
	1.0	M1-2-1410	M1-2-1410-I						
	1.2	M1-2-1412	M1-2-1412-I	M1-2-1812	M1-2-1812-I				
	1.5	M1-2-1415	M1-2-1415-I	M1-2-1815	M1-2-1815-I				
	1.7	M1-2-1417	M1-2-1417-I	M1-2-1817	M1-2-1817-I	M1-2-2017	M1-2-2017-I	M1-2-2417	M1-2-2417-I
	2.0	M1-2-1420	M1-2-1420-I	M1-2-1820	M1-2-1820-I				
	2.4			M1-2-1824	M1-2-1824-I	M1-2-2024	M1-2-2024-I	M1-2-2424	M1-2-2424-I
	2.5	M1-2-1425	M1-2-1425-I						
	3.0	M1-2-1430	M1-2-1430-I	M1-2-1830	M1-2-1830-I				
	3.4			M1-2-1834	M1-2-1834-I	M1-2-2034	M1-2-2034-I	M1-2-2434	M1-2-2434-I
	4.4			M1-2-1844	M1-2-1844-I	M1-2-2044	M1-2-2044-I	M1-2-2444	M1-2-2444-I

Feed Sets Included in Kit



MiniONE® Button comes packaged as qty/box 1 (sterile). Kit contents include: device (1), bolus extension set (1), y-port extension set (1), 10cc syringe (1), 35cc bolus syringe (1), Placement T-Handle/Removal T-Handle[‡], snap arm assembly, and gauze.

Spare Removal Tools For Obturated Buttons

FRENCH SIZE	DESCRIPTION	ORDER #	QTY/BOX
14F	1.0, 1.2 cm	R-S1012	1
14F	1.5, 1.7 cm	R-S1517	1
14F	2.0 cm	R-S2023	1
14F	2.5 cm	R-S2427	1
14F	3.0 cm	R-S3033	1
18F	1.2 cm	R-L1012	1
18F, 20F, 24F	1.5, 1.7 cm	R-L1517	1
18F	2.0 cm	R-L2023	1
18F, 20F, 24F	2.4 cm	R-L2427	1
18F	3.0 cm	R-L3033	1
18F, 20F, 24F	3.4 cm	R-L3437	1
18F, 20F, 24F	4.4 cm	R-L4044	1

AMT Obturator for Distending button during insertion and removal.



[‡]Non-Balloon Kit: Placement and Removal T-Handle
Capsule Non-Balloon Kit: Removal T-Handle

For a complete list of products, visit us at www.AppliedMedical.net
800 869 7382 | CS@AppliedMedical.net | ICS@AppliedMedical.net



Applied Medical Technology, Inc.
8006 Katherine Boulevard
Brecksville, OH 44141
www.AppliedMedical.net
An ISO 13485 Registered Company
Made In USA



© 2026 Applied Medical Technology, Inc. **C-4550-G**
Patents: <http://www.AppliedMedical.net/Patents/> | Trademarks: <http://www.AppliedMedical.net/Trademarks/>
1. Based on internal testing of 14F, 18F, 20F, and 24F devices. | ENFit® is a registered trademark of GEDSA.